



ACROS ORGANICS part of Thermo Fisher Scientific

ACROS
ORGANICS

Product Specifications

Product 19661
Chromium(VI) oxide, 99.5%

General Product Data

Version	4
CAS No	1333-82-0
Molecular weight	99.99
Molecular formula	Cr O ₃
Linear formula	CrO ₃

Product Specifications

Appearance (Color)	Dark red to purple
Appearance (Form)	Crystals or granules
Titration Iodometric	>=99.4 %
Sulfate (SO ₄)	=<500 ppm

Product Specifications are subject to amendment and may change over time.

Acros Organics

ENA23, zone 1, nr 1350, Janssen Pharmaceuticaaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: <http://www.acros.com>
1 Reagent Lane, Fair Lawn, NJ 07410, USA
Fax 201-796-1329