



Certificate of Conformity

VWR FILTER FUNNELS WITH LONG STEM

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	Internal diameter 'D2'	Stem height 'H1'	Stem diameter D1'	Stem thickness 'S'	Pk
221-0265	25±1	150±3	6±0.3	1.5±0.15	10
221-0266	35±1	150±3	8±0.3	1.5±0.15	10
221-0267	50±1	150±3	8±0.3	1.5±0.15	1
221-0268	65±1	150±3	8±0.3	1.5±0.15	10
221-0269	75±1	150±3	8±0.3	1.5±0.15	10
221-0270	100±1.5	150±3	10±0.3	1.5±0.15	1

Material : Borosilicate glass 3.3

Other Characteristics:

- With smooth internal walls.
- Very good chemical resistance.
- High temperature resistance.
- Minimal thermal expansion giving relatively high resistance to temperature changes

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **February 2017**