

TSE/BSE - Declaration

20904 2-PROPANOL PH.EUR./USP

With reference to guideline EMA/410/01 Rev. 3 of the EC, we declare that this product is manufactured without the use of raw materials of animal or human origin.

The product is produced synthetically by catalytic hydration of Propene, followed by distillation.

During processing the product does not come in contact with animal material. Therefore the product is considered safe regarding TSE.



Evelyn Verstraete
QA VWR Chemicals



Esther Venmans
Manager QA
Chemical Manufacturing Europe

19/04/2017