

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Rectangular Micro Cover Glasses
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2021-06-11

No. 1 Slides with thickness 0.13 to 0.17 mm

VWR Catalogue Number	Length x Width [mm]
48393-026	30 x 22
48393-059	50 x 22
48393-070	60 x 22
48393-092	30 x 24
48393-060	40 x 24
48393-081	50 x 24
48393-106	60 x 24

No. 1 - 1/2 Slides with thickness 0.16 to 0.19 mm

VWR Catalogue Number	Length x Width [mm]
48393-172	40 x 22
48393-230	40 x 24
48393-241	50 x 24
48393-251	60 x 24
48404-467	30 x 24
48393-221	60 x 22
48393-151	30 x 22
48393-195	50 x 22

Quality System Compliance

Products are manufactured under the [ISO 13485:2016](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Borosilicate glass
Sterilization	Not Applicable.
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	Not Applicable.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 01/19/2021 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store at general laboratory storage condition
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.