



Certificate of Conformity

VWR MULTI LAYER ADHESIVE FLOOR MATS, PURESTEP

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	L×W (mm)	Thickness	Colour	Layers	Pk
115-0168	660×1143	2.7mm	White	60	4
115-0169	914×1143	2.7mm	White	60	4

Other Specifications:

Film : Low density Polyethylene
Adhesive : Acrylic
Adhesive Thickness : 8 micron
Peel Adhesion : 190g/25mm $\pm$ 40 (individual layers)
Heat Resistance : 70°C - 48 hours

Other Characteristics:

- Bacteriostatic, prevents the growth of microorganisms
- With numbered corner tabs, to know how many layers remain
- Reduced maintenance and cleaning costs

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **April 2017**