

NON STERILE STAINLESS STEEL HANDLES

Supplied in boxes of 10

Handles in this range:

No.3 Graduated, 3L, No.7, No.9
Fit all No.3 fitment blades

No.4 Graduated, 4L
Fit all No.4 fitment blades

COMPONENT	MATERIAL
HANDLE	Stainless Steel 17-4ph
UNIT PACK	Cellophane Sleeve
UNIT CONTAINER BOX	Coated white lined chipboard

Supplied individually in wallets

Handles in this range:

No.5B, B3
Fit all No.3 fitment blades

No.6B
Fit all No.4 fitment blades

COMPONENT	MATERIAL
HANDLE	Stainless Steel 17-4ph
UNIT PACK	PVC Wallet

All handles are manufactured in accordance with BS 2982 and BS EN 27740.

All Components are manufactured in Europe.

All packets and boxes are marked with the current lot number and expiry dates for traceability.