



Model:
CCD114 - 630-1644

Description:
C-Mount adapter for 1/2" sensor, magnification 0,45x

Optical magnification	0,45x
Connection to the camera	C-mount standard
CS-mount	Allowed, but the C-mount ring (optional) is suggested
Connection to the microscope	23mm eyepiece connection
Allowed connections	30mm eyepiece connection (with M113.1)
Sensor size compatibility	1/2" 1/2,3" 1/2,5"
External fast identification	a BLU line
Length	35mm
Color	Black
Allowed	The camera is free to move 360 degree