



## VWR® Traceable® GIANT-DIGIT Countdown Timer

View 3.4 cm giant size digits on your bench timer from across lab

- Features: countdown/up, memory, stopwatch mode
- Unique memory returns channel to previously programmed setting at touch of a key
- Giant-Digits™ with 3.39 cm digit height; Viewable from 15 metres
- Countdown and alarms on single channel for ease-of-use
- Three-button intuitive operation eliminates reading instructions
- Automatic Bounceback™ Memory for repetitive timing
- Counts up during alarming indicating elapsed time since alarming
- Settable to one second with one-second resolution
- Accuracy is 0.01%
- Counts up in stopwatch mode
- Easel bench stand and magnetic back

### Applications

Always on time, timer has a wide variety of ways to suit your laboratory needs.

### Traceable® Certificate

An individually-numbered Traceable® Certificate provided with each unit, assures accuracy from an ISO/IEC 17025:2005 (1750.01) calibration laboratory accredited by A2LA. It indicates traceability of measurements to the SI units through NIST or other recognized national measurement institutes (NMI) that are signatories to the CIPM Mutual Recognition Agreement.

## VWR® Traceable® GIANT-DIGIT Countdown Timer

| Timing Capacity        | Resolution | Accuracy | Cat. No. | UOM   |
|------------------------|------------|----------|----------|-------|
| 99 minutes, 59 seconds | 1 second   | 0.01%    | 609-0231 | 1 ea. |

Individually calibrated ISO 17025 certified

Individually certified—traceable to NIST