

Certificate of Conformity

VWR Catalogue Number	611-4672
Description	VWR® Antistatic Weighing Boat, 85 ml, with 85×85×24 mm Dimension, Medium, Sterile
Lot Number	519259
Date of Manufacture (yyyy-mm-dd)	2021-12
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2023-02-24

Quality System Compliance

Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected in accordance with AQL Level II 4.0 sampling plan per MIL-STD-1916.

Product Specifications

Material	100% Virgin, Vacuum-molded of polystyrene homo-polymer with non-toxic pigment. Materials used meet FDA 21 CFR 177.1640 Compliance.
Temperature Resistance	–10 to 70 °C
Sterilization	Product is Electron Beam irradiated with a specified dose range of 25 – 50 kGy.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	According to EU Directive 2015/863 No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-06-10 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Storage temperatures -20 °C to 60 °C (-4 °F to 140 °F)
Shelf Life	5 Years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.