

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® PHP Plates
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2025-01-02

VWR Catalogue Number	Description	Diameter [mm]	Packed
391-0685	PHP Petri plate, stack size 10	90	10x60
391-0686	PHP Petri plate, stack size 10	90	10x60
391-0687	PHP Petri plate, stack size 10	90	10x60
391-0692	PHP Petri plate, stack size 10	90	10x60
391-0688	PHP Petri plate, stack size 20	90	20x30
391-0689	PHP Petri plate, stack size 20	90	20x30
391-0694	PHP Petri plate, stack size 20	90	20x30

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is manufactured in 10000 cleanroom based on Class 7 specification according to ISO 14644-1 standards.

Material

Polystyrene.
In compliance with 21 CFR 177.1640.

Sterilization	Products are Gamma irradiated and dose released upon ANSI/AAMI/ISO 11137 ISO 13004. Products meets a minimum sterility Assurance Level (SAL) of 10 ⁻⁶ .
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Not Applicable.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	The acceptance level for product is 0.01 EU/ml.(LAL) method
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Material comply with the limit as set by Directive (EU) 205/863 amending Annex II to Directive 2011/65/EU.
REACH Statement	Product complies to >0.1% by weight of the substance of the very high concern as listed by regulation No.EC 1907/2006.
Storage Conditions	Room temperature.
Shelf Life	3 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
 This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	11-2024
2	2.0	ISO, Sterilization, DNase RNase, BSE-TSE, Cytotoxicity, Non-pyrogenic, Latex, REACH, ROHS	01-2025
3	1.2		