

Certificate of Conformity

VWR Catalogue Number	228-0701
Description	Tubing, silicone, 1.00 mm internal diameter, 3.00 mm external diameter, 1.00 mm wall thickness
Lot Number	5324326
Date of Manufacture (yyyy-mm)	2022-06
Expiry Date (yyyy-mm)	2024-12
Country of Origin	Manufactured in France
Date of Issue (yyyy-mm-dd)	2024-02-28

Quality System Compliance

Products are manufactured under the **ISO 13485** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Silicone, Peroxide cross-linked. Formulation: S60V34B. Product complies with the GMP conditions, FDA Regulations CFR 177.2600 USP <88> Class VI requirements and meets EU Pharmacopoeia 3.1.9.
Temperature resistance	- 50 °C to +200 °C (briefly upto +250 °C).
Sterilization	No Information.
ATP Assay	No Information.
DNase & RNase Free	No Information.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin according to the definitions in EMEA/410/01. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted to quality all material resins using ISO 10993-5 standards for cytotoxicity and have been shown to be non-toxic.
Non-Pyrogenic Statement	No Information.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2024-01-23 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Ambient temperature and humidity conditions, protected from sunlight or fluorescent light or UV.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.