

Data sheet

Article: Burette 25 ml class AS, Schellbach
Code no.: 314 01 70

Material:

DURAN (Borosilicat glass 3.3)

Chemical properties:

- Hydrolytic resistance at 98°C: Class 1 acc. ISO 719
- Hydrolytic resistance at 121°C: Class 1 acc. ISO 720
- Acid resistance: acc. ISO 1776
- Alkali resistance: Class 2 acc. ISO 695

Nominal volume:

25 ml

Subdivision:

0,05 ml

Tolerance:

± 0,03 ml

Inner diameter:

7,54 mm

Standard:

ISO 385

Hirschmann Laborgeräte
GmbH & Co.KG



Roland Benz
Quality Manager

Eberstadt, 06.04.2021