

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® Cleanroom Overalls, White Colour
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2021-10-04

Standard overalls with Collar

VWR Catalogue Number	Size
113-1078	M
113-1079	L
113-1080	XL
113-1081	XXL
113-1082	3XL
113-1084	5XL

Overalls with Collar attached non skid boot covers

VWR Catalogue Number	Size
113-1085	M
113-1086	L
113-1087	XL
113-1088	XXL
113-1089	3XL

Hooded overall with attached fluid-impervious boot cover

VWR Catalogue Number	Size
113-9161	S
113-1074	M
113-1075	L
113-1076	XL
113-1077	XXL

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

Bacterial filtration efficiency 99.9%

	Test Method	MD	CD
Break Tensile Properties [N]	ASTM D5034	> 100	> 100
Elongation rate [%]	ASTM D5034	> 50	> 50
Trapezoid tearing strenth(N)	FZ/T01085-1999	>400	>400

Product Specifications

Material	Spunbonded Polypropylene (SF) 63 g/m².
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The above product is latex free.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Room temperature, dry and shady.
Shelf Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.