

Certificate of Conformity

VWR Catalogue Number	See list
Description	VWR® Standard Micro Syringes, 50 mm length, bevelled tip
Country of Origin	Manufactured in Japan
Date of Issue (yyyy-mm-dd)	2023-09-05

Catalogue Number	Description	Volume [µL]
549-0199	Standard micro syringe	10
549-0522	Standard micro syringe	25
549-0523	Standard micro syringe	50
549-0524	Standard micro syringe	100

Quality System Compliance

Products are manufactured under the **ISO 14001:9001** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Borosilicate glass, Stainless steel, Ni coated brass
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Not Applicable.
Cytotoxicity	Not Applicable.

Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2023-06-14 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Temperature -10°C to 50°C. Humidity 20% to 90%
Shelf Life	5 years since manufacture. 1 year after delivery.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	09-2023