

Certificate of Conformity

VWR Catalogue Number	See List
Description	Snap Cap Vial
Country of Origin	Manufactured in India - 548-0141A Manufactured in China - 548-0142A, 548-0143A, 548-0554A, 548-0555A, 548-0625A, 548-0620A, 548-0621A, 548-0623A, 548-0649A, 548-0622A, 548-0624A, 548-0650A, 548-0651A
Date of Issue (yyyy-mm-dd)	2022-09-14

Vials, Flat Bottom, Borosilicate glass 5.1, Clear, without Cap

Catalogue Number	Capacity [mL]	Dimension [mm]	Cap Size
548-0141A	5.00	20 × 40	ND18
548-0142A	10.00	22 × 50	ND18
548-0143A	15.00	26 × 48	ND22

Vial Kits, Flat bottom, Clear, Polyethylene Cap, Transparent Cap Colour

Catalogue Number	Capacity [mL]	Dimension [mm]	Cap Size
548-0554A	3.00	30 × 18	ND18
548-0555A	5.00	20 × 40	ND18
548-0625A	11.00	22 × 45	ND22
548-0620A	15.00	24 × 52	ND22
548-0621A	12.50	22 × 50	ND22
548-0623A	20.00	26 × 55	ND22
548-0649A	25.00	30 × 50	ND28
548-0622A	27.00	26 × 70	ND22
548-0624A	30.00	28 × 75	ND22
548-0650A	40.00	30 × 80	ND28

Catalogue Number	Capacity [mL]	Dimension [mm]	Cap Size
548-0651A	50.00	30 × 100	ND28
548-0652A	64.00	28 × 145	ND28

Quality System Compliance

Products are manufactured under the [ISO 9001:2015](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Borosilicate glass 5.1. Polyethylene cap.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.

REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Not applicable.
Shelf Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature