

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® PHP Plates
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2024-11-18

VWR Catalogue Number	Description	Diameter [mm]	Packed
391-0685	PHP Petri plate, stack size 10	90	10×60
391-0686	PHP Petri plate, stack size 10	90	10×60
391-0687	PHP Petri plate, stack size 10	90	10×60
391-0692	PHP Petri plate, stack size 10	90	10×60
391-0688	PHP Petri plate, stack size 20	90	20×30
391-0689	PHP Petri plate, stack size 20	90	20×30
391-0694	PHP Petri plate, stack size 20	90	20×30

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material Virgin Crystal Polystyrene; US FDA approved

Dimensional Specification

Base External Diameter	87 ± 1 mm
Base Height	13.70 ± 1 mm
Lid External Diameter	91.70 ± 1 mm
Lid Height	7.52 ± 1 mm
Assembled Height	14.65 ± 1 mm
Assembled Weight	14.75 ± 1 g

VWR Singapore Pte Ltd, 18 Gul Drive, Singapore 629468 - <http://www.vwr.com>

Sterilization	Products are Gamma irradiated as per ISO 11137 which meets a minimum requirement of 10 ⁻⁶ SAL with a specified dose range of 9-15 kGy
ATP Assay	Not Applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	This product is free of cytotoxicity effects.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Room temperature.
Shelf Life	3 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	11-2024
2	1.1		
3	1.2		