

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Centrifuge Tubes, ASTM
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2025-04-28

Catalogue Number	Description	Capacity [mL]	Ø×L [mm]	Closure Type	Pk
525-1568	Centrifuge tube, conical	100	17 × 200	Without cap	1
525-1569	Centrifuge tube, pear-shaped	100	17 × 200	Without cap	1

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 17025:2017** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is in conformance with the requirements of the following standards and regulations

Technical standards for products: ASTM E-982.

Calibration procedure: NA

Material The products are made from Borosilicate Glass 3.3 Expansion.

Sterilization Hot air sterilization, in the oven up to 140 °C
Steam sterilization, in an autoclave 121 °C/ 20 min
134 °C/ 10 min.

ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable
Storage Conditions	Store at room temperature in normal warehouse conditions.
Shelf Life	No predefined shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	04-2025
2	1.1	Modified Sterilization information	04-2025