

**Certificate of Analysis**

**Product: STERILE WATER FOR IRRIGATION, USP - 1000 ML**

**Lot #: G116467**  
**Code: 2F7114**  
**Expiry: Feb 2018**

**Manufactured Date: 02/13/2015**  
**Sterilization Date: 02/13/2015**  
**Mixing Specification: 03-15-53-216**  
**Issue Date 11/13/2007**

| TEST NUMBER  | TEST NAME                                                   | LIMIT                            | RESULT  |
|--------------|-------------------------------------------------------------|----------------------------------|---------|
| 11-20-15-001 | Oxidizable Substances                                       | Pass                             | Pass    |
| 11-21-16-003 | pH at 25 deg. C                                             | 5.0 - 7.0                        | 5.6     |
| 11-22-03-007 | Conductivity Determination of Water (USP/EP)                | NMT 5 uS/cm                      | 1 uS/cm |
| 20-03-07-006 | Particle Analysis - Large Volume (Instrumental)             |                                  |         |
|              | Particles >= 10 micron                                      | NMT 25 particles/mL >=10 microns | 0 ct/mL |
|              | Particles >= 25 micron                                      | NMT 3 particles/mL >=25 microns  | 0 ct/mL |
| 13-01-V      | Bacterial Endotoxin                                         | < 0.25 EU/mL = Pass              | Pass    |
| 14-04-01-120 | Sterility                                                   | Pass Parametric Release          | Pass    |
|              | PARAMETRIC RELEASE IN ACCORDANCE WITH U.S. FDA REQUIREMENTS |                                  |         |

THIS PRODUCT WAS MANUFACTURED ACCORDING TO CURRENT GMP,  
FULLFILLS THE REQUIREMENTS OF THE MASTER PRODUCTION DOCUMENT AND TESTED TO USP COMPLIANCE.

*S. Kitchens 2/17/15*

Prepared By/Date

SHARON KITCHENS

*Twila Proctor 2/18/15*

Quality Management Approval /Date

Twila Proctor