

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® Screw Thread Headspace Vials, ND18
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2024-10-08

Precision screw headspace vial, Borosilicate glass, No cap included, ND 18 Cap size

Catalogue Number	Capacity [mL]	Ø × H [mm]	Colour	Cap type	Pk
VWRA548-2678	20	22.5 x 75.5	Clear	Screw	600
VWRA548-2679	20	22.5 x 75.5	Amber	Screw	500

Screw thread caps with septa

Catalogue Number	Description	Type	Septum	Ø [mm]	Pk
VWRA548-2680	Screw polypropylene metal cap	Central hole 8 mm	Blue PTFE/white silicone	18	5.000
VWRA548-2681	Screw polypropylene metal cap	Central hole 8 mm	Gray PTFE/Red Butyl	18	5.000
VWRA548-2682	Magnetic precision screw metal cap	Central hole 8 mm	Red PTFE/white silicone	18	5.000
VWRA548-2683	Magnetic precision screw metal cap	Central hole 8 mm	White PTFE/blue transparent silicone	18	5.000

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & GB/T 19001-2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Vials: Borosilicate glass Cap: Metal cap Septa: Refer table.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE) according to EMEA 410/01 revision 3.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Raw materials comply with the requirements of the Directive 2002/95/EC, as amended, and 2011/65/Eu concerning the limits of cadmium, lead, mercury, hexavalent chromium, Polybrominated Biphenyls (PBB) and Polybrominated Diphenyl Ethers (PBDE).
REACH Statement	Not Applicable.
Storage Conditions	Store at room temperature, away from direct sunlight and rain.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	10-2024