


VWR® Slide Zipper bags



Description: VWR® Slider Zipper bag, LDPE, transparent, with blue slider

Purpose: Packing goods

Materials: LDPE 260GG+FD21SH

Features:

- Uses: under -10 ~65 °C
- Sterilization: N/A
- BSE/TSE: The above product contains no ingredients of animal origin and no material derived from, or exposed to animals affected by or under quarantine for Transmitting, Animal Spongiform Encephalopathy (TSE)/Bovine Spongiform Encephalopathy (BSE).
- Products are manufactured under the ISO 9001:2015 standard.
- Shelf Life: N/A.

VWR® Slider Zipper Bag

EU Part No.	Product Description		
VWRI129-1200	VWR® Slider Zipper bag, LDPE, transparent,250mmx350mmx50um, Write on, with blue slider		
VWRI129-1201	VWR® Slider Zipper bag, LDPE, transparent,300mmx400mmx50um, Write on, with blue slider		

EU Part No.	Width(mm)	Length(mm)	Thickness(um)	Printing	Printing size	Side strength	Tensile strength		Elongation
	Tolerance(+3%/-2%)	Tolerance(+3%/-2%)	Tolerance(±5%)				MD	TD	
VWRI129-1200	250	350	50	Yes	210mmx250mm	10N/15mm	21	21	550%
VWRI129-1201	300	400	50	Yes	240mmx300mm	10N/15mm	21	21	550%