




Europe Cat#	NA Cat#	Description	Material	Capacity	Graduation	Product Dia (mm)	Product Max dimension with Spout & handle (mm)	Product Height (mm)	Net Weight (Grams)
213-0497	75800-250	VWR® Beaker With Handle 100ml - PP	Polypropylene	100ml	10 ml	63.40	95.00	74.00	23.00
213-0498	75800-252	VWR® Beaker With Handle 250ml - PP	Polypropylene	250ml	25 ml	84.00	116.00	94.00	48.00
213-0499	75800-254	VWR® Beaker With Handle 500ml - PP	Polypropylene	500ml	50 ml	110.00	152.00	128.50	88.00
213-0500	75800-256	VWR® Beaker With Handle 1000ml - PP	Polypropylene	1000ml	50 ml	129.00	175.00	151.00	161.00
213-0501	75800-258	VWR® Beaker With Handle 2000ml - PP	Polypropylene	2000ml	100 ml	158.00	216.00	192.30	245.20
213-0502	75800-260	VWR® Beaker With Handle 3000ml - PP	Polypropylene	3000ml	200 ml	177.00	220.00	216.00	305.50
213-0503	75800-262	VWR® Beaker With Handle 5000ml - PP	Polypropylene	5000ml	500 ml	203.00	260.00	262.70	423.20
213-0504	76019-278	VWR® Beaker With Handle 100ml - TPX	TPX / PMP	100ml	10 ml	63.40	96.00	73.50	22.50
213-0505	76019-280	VWR® Beaker With Handle 250ml - TPX	TPX / PMP	250ml	25 ml	84.00	116.50	93.00	47.00
213-0506	76019-282	VWR® Beaker With Handle 500ml - TPX	TPX / PMP	500ml	50 ml	110.50	153.00	128.00	87.00
213-0507	76019-284	VWR® Beaker With Handle 1000ml - TPX	TPX / PMP	1000ml	50 ml	129.40	176.00	150.00	160.50
213-0508	76019-286	VWR® Beaker With Handle 2000ml - TPX	TPX / PMP	2000ml	100 ml	158.00	217.00	191.00	245.00
213-0509	76019-288	VWR® Beaker With Handle 5000ml - TPX	TPX / PMP	5000ml	500 ml	203.00	261.00	262.00	422.00



Material Recycle Code



PP



OTHER

Caution:

1. Do not stack products while Autoclaving, this may lead to physical damage / collapse of product.

Material: Refer above table for material type, 100 % Virgin, FDA 21 CFR Compliance.

Description: Beaker With handle

Drawing:
Not to Scale



Release: 08/30/17

Revision: 0