

VWR® Spec-Wipe 4 Wiper

Made of 100% continuous filament polyester, delivering softness, high absorbency, strength, and abrasion resistance. The wipes are laundered in an ISO Class 4 (FED-STD-209E Class 10/M2.5) cleanroom. The wipes are packaged in two layers of cleanroom grade, low-density polyethylene bags. The inner bags are labeled with a plain black and white label. In Addition, wipes are available sterile by Sterility Assurance Level (SAL) of 10^{-6} in accordance with ANSI/AAMI/ISO 11137 guidelines.

- Continuous filament
- Textured yarn construction
- Clean-cut edges

Technical Data

VWR® Spec-Wipe® 4 Wiper		
Attribute; (units)	Value	Test Method
Basis weight; (g/m ²)	136	Contec Method*
Sorbency in water		
Extrinsic; (mL/m ²)	312	IENT-RP-CC004.3, Sec. 8.1
Sorptive rate; (seconds)	<1.5	
Extractables		IENT-RP-CC004.3, Sec. 7.1.2
In deionized water; (g/m ²)	0.042	
In isopropyl alcohol; (g/m ²)	0.015	
Specific ions		IENT-RP-CC004.3, Sec. 7.2.2
Sodium; (ppm)	0.39	
Chloride (ppm)	0.45	
Particle generation		IENT-RP-CC004.2, Sec. 5.1
P ₀ ≥ 0.5µm; (x10 ⁶ /m ²)	4.8	
Fibers > 100µm; (x10 ³ /m ²)	3.3	IENT-RP-CC004.2, Sec. 5.2

Part Number	Description	EA/BG	BG/CS
76377-480	Spec-Wipe 4 CS sterile continuous filament polyester, 30.5 x 30.5 cm (12 x 12")	300	1