




Europe Cat#	NA Cat#	Description	Material	Capacity	Graduation	Product Dia (mm)	Product Max dimension with Spout (mm)	Product Height (mm)	Net Weight (Grams)
213-0510	75800-234	VWR® Beaker without handle, 25ml - TPX	TPX / PMP	25ml	5 ml	43.0	46.3	50.0	7.1
213-0482	75800-236	VWR® Beaker without handle, 50ml - TPX	TPX / PMP	50ml	5 ml	50.5	54.0	58.0	10.2
213-0483	75800-238	VWR® Beaker without handle, 100ml - TPX	TPX / PMP	100ml	10 ml	63.4	71.0	73.5	18.0
213-0484	75800-240	VWR® Beaker without handle, 250ml - TPX	TPX / PMP	250ml	25 ml	84.0	91.0	93.0	42.1
213-0485	75800-242	VWR® Beaker without handle, 500ml - TPX	TPX / PMP	500ml	50 ml	110.5	118.0	128.0	77.0
213-0486	75800-244	VWR® Beaker without handle, 1000ml - TPX	TPX / PMP	1000ml	50 ml	129.4	140.6	150.0	131.0
213-0487	75800-246	VWR® Beaker without handle, 2000ml - TPX	TPX / PMP	2000ml	100 ml	158.0	172.0	192.0	202.0
213-0488	75800-248	VWR® Beaker without handle, 5000ml - TPX	TPX / PMP	5000ml	500 ml	203.0	223.0	262.0	345.0
213-0511	76019-260	VWR® Beaker with out Handle 25ml - PP	Polypropylene	25ml	5 ml	42.5	46.5	49.8	7.4
213-0489	76019-262	VWR® Beaker with out Handle 50ml - PP	Polypropylene	50ml	5 ml	49.8	54.0	57.7	11.0
213-0490	76019-264	VWR® Beaker with out Handle 100ml - PP	Polypropylene	100ml	10 ml	63.4	70.0	72.5	19.0
213-0491	76019-266	VWR® Beaker with out Handle 250ml - PP	Polypropylene	250ml	25 ml	84.0	90.5	92.5	44.0
213-0492	76019-268	VWR® Beaker with out Handle 500ml - PP	Polypropylene	500ml	50 ml	110.0	117.0	127.5	81.0
213-0493	76019-270	VWR® Beaker with out Handle 1000ml - PP	Polypropylene	1000ml	50 ml	129.0	140.0	149.5	139.0
213-0494	76019-272	VWR® Beaker with out Handle 2000ml - PP	Polypropylene	2000ml	100 ml	158.0	171.0	191.0	215.0
213-0495	76019-274	VWR® Beaker with out Handle 3000ml - PP	Polypropylene	3000ml	200 ml	177.0	193.0	230.0	292.0
213-0496	76019-276	VWR® Beaker with out Handle 5000ml - PP	Polypropylene	5000ml	500 ml	203.0	222.0	261.0	361.0

Material Recycle Code



Caution:

1. Do not stack products while Autoclaving, this may lead to physical damage / collapse of product.

Material: Refer above table for material type, 100 % Virgin, FDA 21 CFR Compliance.

Description: Beaker Without Handle

Drawing:
Not to Scale



Release: 08/30/17

Revision: 0