

Certificate of Conformity

VWR Catalogue Number	See List
Description	Overalls, with hood, VWR® Basic
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2022-03-02

Catalogue Number	Clothing Size	Colour
113-9026	S	White
113-9027	M	White
113-9028	L	White
113-9029	XL	White
113-9030	XXL	White
113-9031	3XL	White
113-9032	S	Blue
113-9033	M	Blue
113-9034	L	Blue
113-9035	XL	Blue
113-9036	XXL	Blue
113-9037	3XL	Blue

Quality System Compliance

Products are manufactured under the **ISO 13485:2012** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

	Test Method	MD	CD
Break Tensile Properties [N]	GB/T3923.1-1997	> 50	> 40
Elongation rate [%]	GB/T3923.1-1997	> 30	> 30

Product Specifications

Material	40 gsm Polypropylene (PP) laminated with Polyethylene (PE) film.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Room temperature, dry and shady.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

VWR International bv, Haasrode Researchpark 2020, Geldenaaksebaan 464, B-3001 Leuven, Belgium - <http://www.vwr.com>