

Precellys Emulsion Kit

Get Perfect Emulsions



- **Easy:** put your Adjuvant and Antigen (peptide, bacteria...) solutions in the tube and let the instrument mix for you
- **Fast & Efficient:** only 2 min shaking time on Minilys for a perfect emulsion
- **Flexible:** prepare 2 to 8mL of emulsion
- **Safe:** no risk of accidental exposure thanks to the single tube concept and no Repetitive Stress Injuries (RSI)

An efficient method...



Precellys Emulsion Kit




PERFECT EMULSIONS




...to prepare homogenous and stable emulsions


Precellys
Emulsion kit



2-syringe method



Vortex



Precellys Emulsion Kit can be used for:

- Preparation of emulsions with Freund's Adjuvant + Antigens
- Monoclonal & polyclonal antibody production
- Preclinical animal models for research & in industry
- DTH experiments (Delayed Time Hypersensitivity)

REFERENCE: D34200

Content

- 10 autoclaved pouches including:
- 1x Emulsion tube and cap
 - 1x specific Emulsion tube plunger

DISCOVER OUR LIFE SCIENCES PRODUCT RANGE



AIR SAMPLING



SAMPLE
PREPARATION



ASSAY KITS



CELLIMAGING
SYSTEM

Bertin Technologies - HeadOffice

10 bis Avenue Ampère - Parc d'activités du Pas du Lac
78180 Montigny-le-Bretonneux FRANCE