

Terumo Syringe

3-part Syringe



Product Information

NEW PRODUCT CODES!



Disposable syringes for single use

Every single detail is the result of Terumo's long time experience in this field which is appreciated by users throughout the world.

- Siliconisation for optimal plunger mobility
- The integral gasket ensures a firm and distinct end-stop
- Clear graduated scale allows accurate dosage
- Transparent barrel for visualisation
- Finger grips to support easy handling

Description	Raw Materials
1 Barrel	Polypropylene
2 Plunger	Polypropylene
3 Gasket	Thermoplastic elastomer (for ,SS'-codes) Polyisoprene rubber (for ,MDSS'-codes)
Lubricant	Silicone oil

Does not contain components made of PVC or natural rubber latex.



Sterilization

With radiation (e-beam) (for ,SS'-codes) or ethylene oxide (for ,MDSS'-codes) to a Sterility Assurance Level of at least 10^{-6}



Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the paper (for ,MDSS'-codes) or embossed in the film of the blister pack (for ,SS'-codes)
- Do not store at extreme temperature and humidity



Precautions

- Do not use if package is damaged
- Do not re-use
- Dispose of sharps after single use according to facility protocol to avoid risk of infection

Always consult the product label for a complete overview of warnings, cautions and/or precautions prior to actual use.

 **TERUMO**

Your syringes of choice in different sizes and tip types


NEW PRODUCT CODES!

Product Range

Product Code	Nominal Capacity	Tip	Replacement for
MDSS01SE	1 ml	Concentric Luer Slip	SS+01T1
SS*02LE1	2,5 ml	Concentric Luer Lock	-
SS*02SE1	2,5 ml	Concentric Luer Slip	-
MDSS03SE	3 ml	Concentric Luer Slip	SS*02SE1, SS+03S
MDSS03LE	3 ml	Concentric Luer Lock	SS*03L, SS+03L1
MDSS05SE	5 ml	Concentric Luer Slip	SS*05SE1, SS+05S
SS*05LE1	5 ml	Concentric Luer Lock	-
SS*05SE1	5 ml	Concentric Luer Slip	-
MDSS10SE	10 ml	Concentric Luer Slip	SS+10ES1
MDSS10ESE	10 ml	Eccentric Luer Slip	SS+10ES1
SS*10LE1	10 ml	Concentric Luer Lock	-
MDSS20ESE	20 ml	Eccentric Luer Slip	SS+20ES1
SS+20L1	20 ml	Concentric Luer Lock	-
SS+30L1MP	30 ml	Concentric Luer Lock	-
SS+30ES1	30 ml	Eccentric Luer Slip	-
SS+50L1	50 ml	Concentric Luer Lock	-
SS+50ES1	50 ml	Eccentric Luer Slip	-
SS+50C1	50 ml	Catheter Tip	-

Packaging

Unit Box	Shipping Carton		
	Quantity	Gross Weight	Dimensions (L x W x H)
100	3000	16.60 kg	635 x 295 x 505 mm
100	2400	6.00 kg	397 x 581 x 444 mm
100	2400	6.00 kg	394 x 581 x 444 mm
100	1800	10.90 kg	445 x 447 x 449 mm
100	1800	10.90 kg	445 x 447 x 449 mm
100	1800	14.90 kg	656 x 419 x 426 mm
100	1600	6.25 kg	412 x 573 x 444 mm
100	1600	6.25 kg	411 x 576 x 444 mm
100	1200	14.60 kg	569 x 473 x 459 mm
100	1200	14.60 kg	565 x 425 x 485 mm
100	800	15.00 kg	359 x 581 x 383 mm
50	400	14.60 kg	468 x 390 x 455 mm
50	600	10.15 kg	573 x 410 x 377 mm
25	200	6.32 kg	560 x 398 x 255 mm
25	200	6.20 kg	558 x 395 x 254 mm
25	100	5.10 kg	464 x 320 x 275 mm
25	100	5.10 kg	472 x 324 x 282 mm
25	100	5.40 kg	472 x 324 x 282 mm



Single unit: Blister pack
(film/film for ,SS'-codes and film/paper for ,MDSS'-codes)

Regulatory Status

CE Certification – Class IIa medical device

Notified body TÜV Rheinland 0197

Manufacturer for ,SS+/-codes TERUMO (Philippines) CORPORATION,
124 East Main Avenue, Laguna Technopark, Binan, Philippines
for ,SS*'-codes TERUMO CORPORATION, 2-44-1 Hatagaya,
Shibuya-Ku, Toyko 151-0072, Japan
for ,MDSS'-codes Yangzhou Medline Industry Co., Ltd. (China)

EC representative for ,SS'-codes TERUMO EUROPE N.V. – Leuven, Belgium
for ,MDSS'-codes Shanghai International Corp. GmbH (Europe) – Hamburg, Germany

Distributor TERUMO EUROPE N.V., Interleuvenlaan 40, 3001 Leuven, Belgium

Choose from a
range of volumes
from 1 ml to 50 ml

Syringes are available in all
common tip types.



Terumo Corporation
+81 3 6742 8500

Terumo Europe N.V.
+32 16 38 12 11

Your Terumo partner

All brand names are trademarks
or registered trademarks of
TERUMO CORPORATION and their
respective owners.
Published by Terumo Europe N.V.

Please note that some products may not
be available in all regions, check with
your local Terumo contact.