

Technical Datasheet

VWR WIPES BLUEPROP

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Applicable part numbers : 115-2418, 115-2419

Attribute	Value*			Test method
Composition :				
polypropylene	100%			
Production method	meltblown			
Colour	blue			
Basis weight	70 g/m ²			
Absorption :				
water	660%			
oil	1350%			
Tensile strength :	Dry	Wet	Acetone	DIN/EN 29073 T3
Machine direction	30,2 N	30,9 N	29,4 N	
Cross direction	26,2 N	27,9 N	24,1 N	

Characteristics :

Absorbs water and oil ; Resistant to acids and solvents ; Soft ; Synthetic ; Colorfast.
Recyclable. Food Contact Certified. Silicone free.

Applications :

Aerospace, Automotive, High-Tech, Food & Beverage, Laboratories, Printing, Janitorial.

Formats & Packaging :

Cat No	Type	Format size	Quantity	Primary pkg	Secondary pkg	Pallet
115-2418	Industrial roll	38x32cm	1 roll x 50 pcs	Plastic film	n/a	40 rolls
115-2419	Z-folded wipes	38x40cm	20 packs x 25 pcs	Plastic bag	Cardboard case	28 cases

*** Disclaimer:** The information presented here is the result of testing and other evaluations conducted during the development of the product. This information does not constitute a specification.

Date of issue: **March 2017**