

Certificate of Conformity Finntip™ Sterile

Product Code: 9401113
Product Name: Finntip 1000 sterile
Lot Number: 7260U0
Expiry Date: 2022-09

We hereby confirm that these high quality Finntip™ pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of the following contaminants:

Human DNA	$< 1.9 \times 10^{-11}$ g/tip
DNase	$< 9.4 \times 10^{-4}$ U/tip
RNase	$< 6.5 \times 10^{-8}$ U/tip
Endotoxin	< 0.05 EU/ml

The sterility (SAL 10⁻⁶) is ensured by irradiating this product according to ISO 11137 guidelines.

In Joensuu, 4th January 2018



Heidi Salokangas, Quality Manager
Liquid Handling Consumables
Thermo Fisher Scientific Oy