



VWR® 48 Well, No Lid, Non-treated Tissue Culture Plates



Description: 48 Well Tissue Culture Plate, No Lid, Non-treated, Non-Sterilized

Materials

Plate: GPPS (General polystyrene) Colour: Clear

Features

- Temperature: -20°C~+50°C
- Non-pyrogenic
- Available in bulk PE bag
- Autoclavability: No
- Max. volume (single well): 1.55 ml
- Working volume (single well): 0.19~0.29 ml
- Growth area (single well): 0.83 cm²
- Every product is printed with a batch number for traceability
- Manufactured in a class 100,000 room environment
- Manufactured under EN ISO13485 :2003 and ISO9001:2008 quality management system



VWR® 48 Well, No Lid, Non-treated Tissue Culture Plates

VWR EU Cat. No.	Description	Sterile	Qty/Pk	Pk/Box	Box/Case	Qty/Case
734-2888	48 Well Tissue Culture Plate, No Lid, Non-treated, Non-Sterilized	No	50	1	4	200

Cautionary Stements

Ximus. Ximporem fugiaerae vel is dolluptaepta volupti dolore, undam, quatur, ationsequo que aliqui ne pe exeribust, omnimpel ipsamus que labo. Lo et mint estrum, of- ficatio te doluptate coresciis moloritio que consequo quis sit atus, eume nistem quas moloria qui si aut qui tem estis aut liasimiligni