

Certificate of Conformity

Product Information

Product Code: 10347513

Product Name: Sharkskin 150mm 100pk

Lot Number: 116501

QC Data/Specification/Result

Customer	MH011680	GLOBAL LIFE SCIENCES	Grade	MHPM0992	ProdName	15 CM DIA Shark Skin (TM)						
Product	100/PK-MADE IN USA	C.ProdCode	10347513		Width	965.20	mm	Weight	49.14 g/m2	Color	WHITE	
	MD KLEMM RAPIDITY	THICK CD WET T	DENSO	WEIGHT	THICK	MULLEN	FORM	WET BUR	WEIGHT			
	1/16'S I	MLS/MIN	MIL	LBS/11INC	SECONDS	GSM	MICRONS	PSI	NUMBER (	INCHES O	LBS 1389	
116501/003	18.0	130.0	6.4	2.4	22.2	48.0	162.0	20.0	7.0	84.5	13.7	

Customer	MH011680	GLOBAL LIFE SCIENCES	Grade	MHPM0992	ProdName	38.5 CM DIA Shark Skin (TM)						
Product	100/PK-MADE IN USA	C.ProdCode	10347523		Width	965.20	mm	Weight	49.14 g/m2	Color	WHITE	
	MD KLEMM RAPIDITY	THICK CD WET T	DENSO	WEIGHT	THICK	MULLEN	FORM	WET BUR	WEIGHT			
	1/16'S I	MLS/MIN	MIL	LBS/11INC	SECONDS	GSM	MICRONS	PSI	NUMBER (	INCHES O	LBS 1389	
116111/003	16.0	120.0	6.6	2.5	21.1	48.0	167.7	15.0	7.0	102.0	13.7	
116501/003	18.0	130.0	6.4	2.4	22.2	48.0	162.0	20.0	7.0	84.5	13.7	

Customer	MH011680	GLOBAL LIFE SCIENCES	Grade	MHPM0992	ProdName	34 CM DIA SHARKSKIN						
Product	100/PK-MADE IN USA	C.ProdCode			Width	965.20	mm	Weight	49.14 g/m2	Color	WHITE	
	RAPIDITYMD KLEMM CD WET T	THICK	DENSO	WEIGHT	WEIGHT	FORM	MULLEN	WET BUR	THICK			
	MLS/MIN	1/16'S I	LBS/11INC	MIL	SECONDS	GSM	LBS 1389	NUMBER (	PSI	INCHES O	MICRONS	
116111/004	131.7	16.0	1.6	6.4	21.4	51.0	14.5	7.0	20.0	89.5	161.3	
116501/003	130.0	18.0	2.4	6.4	22.2	48.0	13.7	7.0	20.0	84.5	162.0	

Conformance & Quality systems statement

This is to certify that this product conforms to Cytiva specifications.

All products are manufactured via a Quality System certified to **BS EN ISO9001:2015** and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Electronic signature

This document has been electronically produced and is valid without a signature.