

Certificate of Conformity

VWR Magnetic Stirring Bar, Micro, PTFE

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Lot No.	Length(mm)	Ø (mm)	Pk
442-0368	J104403	20	3	10

Date of Manufacture: February 2021

Material: PTFE Covered, Strong Alnico V magnetic core

Properties: Good chemical resistance

Temperature Resistance: -200 °C up to +280 °C

3M® DYNEON® TF1620 – Compliance Status

FDA – Directive FDA 21 CFR 177.1550

USP Class VI – USP 27, NF 22, 2004 for Plastics at 70°C.

RoHS 2 – Directive 2011/65/EU & Dangerous Substances – 76/769/EC

RoHS 2 – Phthalates – Directive 2011/65/EU as amended by 2015/863/EU

EU REACH Regulations – SVHC Current Candidate List

Bisphenol A – Absent

Latex – Absent

ADCF – Animal Derived Components Free (BSE/TSE)- EMEA/410/01/rev.3. GMO Free

Endotoxins – Free

POP, PFOS, PFOA including Homologues, Derivatives and Salts of APFO, DEHIP, Phthalates, VOC Free

ODC – Class 1

Dodd-Frank Act – Gold, Coltan, Cassiterite, Wolframite and their derivatives originating from the Democratic Republic of Congo - Absent

California Proposition 65 – Absent of chemicals contained on the California Proposition 65 list.

It is the responsibility of the user to ensure that the products are suitable for their intended use.

This document has been produced electronically and is valid without a signature

Date of issue: **March 2021**