

VELP Scientifica Solutions for Jar/Leaching Test



JLT 4



JLT 6



FC4S



FC6S



FP4

Number of positions	4	6	4	6	4 (portable)
Independent positions			•	•	
Speed settings rpm	from 10 to 300	from 10 to 300	10-15-30-45-60-90-120-150-200-300	10-15-30-45-60-90-120-150-200-300	from 10 to 300
Time settings	1÷999 min or 1÷99 hours continuous	1÷999 min or 1÷99 hours continuous	- continuous	- continuous	1min – 99h 59min continuous
Lighting system	illuminated back panel	illuminated back panel	illuminated back panel	illuminated back panel	central light
Dimensions (WxHxD) mm	655x404x296 mm	935x404x296 mm	655x404x296 mm	935x404x296 mm	250x310x265 mm
(in)	(25.7x15.9x11.6 in)	(36.8x15.9x11.6 in)	(25.7x15.9x11.6 in)	(36.8x15.9x11.6 in)	(9.8x12.2x10.4 in)
Weight					
Kg	13.0	17.0	12.5	18.0	3.8 kg
(lb)	(28.6)	(37.4)	(27.5)	(39.6)	(8.3 lb)
Power Supply	100÷240 V	100÷240 V	100÷240 V	100÷240 V	100÷240 V
Power	19 W	19 W	26 W	32 W	6 W



www.velp.com