

VWR[®] PUREImage cleanroom paper and notebooks



PUREImage

The most advanced line of clean paper and documentation available for the most demanding applications.



PUREIMAGE CLEANROOM PAPER


For writing, printing and photocopying in cleanroom environments.

Characteristics:

- Dual binding technology
- Suitable for use in class 10 to 1000 (ISO class 4 to 6) cleanrooms
- PUREImage premier paper produced in class 10 (ISO 4) cleanroom
- Excellent printability and heat toner adhesion
- Double sided web cleaned prior to sheeting
- High opacity and durability

Specifications	80 gsm	105 gsm
Basis weight (g/m²)	78–83	102–107
Caliper (µm)	94–104	110–120
Tear strength MD/CD	74/75	80/81
Brightness (%)	88–92	88–92

Description	Colour	Weight	Pk	Cat. No.
Non latex paper, A4 size	White	80 g/m²	2,500	115-0149
Non latex paper, A4 size	Blue	80 g/m²	2,500	115-0150
Non latex paper, A4 size	Violet	80 g/m²	2,500	115-0155
Non latex paper, A4 size	White	105 g/m²	2,500	115-0156
Latex-impregnated paper, A2 size	White	104 g/m²	250	115-1961
Latex-impregnated paper, A3 size	Blue	85 g/m²	1,250	115-1952
Latex-impregnated paper, A3 size	White	85 g/m²	1,250	115-1958
Latex-impregnated paper, A3 size	Blue	104 g/m²	1,250	115-1954
Latex-impregnated paper, A3 size	White	104 g/m²	1,250	115-1959
Latex-impregnated paper, A4 size	White	80 g/m²	2,500	115-0163
Latex-impregnated paper, A4 size	Blue	80 g/m²	2,500	115-0164
Latex-impregnated paper, A4 size	Yellow	80 g/m²	2,500	115-0165
Latex-impregnated paper, A4 size	Green	80 g/m²	2,500	115-0166
Latex-impregnated paper, A4 size	Pink	85 g/m²	2,500	115-1955
Latex-impregnated paper, A4 size	White	104 g/m²	2,500	115-1956
Latex-impregnated paper, A4 size	White	104 g/m²	2,500	115-1957
Latex-impregnated paper, A4 size	Blue	105 g/m²	2,500	115-1951
Latex-impregnated paper, A5 size	White	104 g/m²	5,000	115-1960
Latex-impregnated paper, A5 size	Blue	105 g/m²	5,000	115-1953



CRITICAL PRINT CLEANROOM NOTEBOOK

For writing in cleanroom environments

Characteristics:

- Cleanroom processed and packed
- Amide- and silicone-free
- PE impregnated paper without any organic fillers
- Recyclable
- Custom printing available
- Ruled
- See-through cover

Specifications	105 gsm
Basis weight (g/m²)	102–107
Caliper (µm)	110–120
Tear strength MD/CD	80/81
Brightness (%)	88–92

Size	Colour	No. of pages	pk	Cat. No.
A4	White	100	10	818-0126
A5	White	100	20	818-0127
A6	White	100	40	818-0128

Setting science in motion to create a better world



AUSTRIA

VWR International GmbH
Graumannsgasse 7
1150 Wien
Tel.: +43 1 97 002 0
info.at@vwr.com

BELGIUM

VWR International bv
Researchpark Haasrode 2020
Geldenaaksebaan 464
3001 Leuven
Tel.: +32 (0) 16 385 011
vwr.be@vwr.com

CZECH REPUBLIC

VWR International s. r. o.
Veetee Business Park
Pražská 442
CZ - 281 67 Stříbrná Skalice
Tel.: +420 321 570 321
info.cz@vwr.com

DENMARK

VWR International A/S
Tobaksvejen 21
2860 Søborg
Tel.: +45 43 86 87 88
info.dk@vwr.com

FINLAND

VWR International Oy
Valimotie 17-19
00380 Helsinki
Tel.: +358 (0) 9 80 45 51
info.fi@vwr.com

FRANCE

VWR International S.A.S.
Immeuble Estréo
1-3 rue d'Aurion
93114 Rosny-sous-Bois cedex
Tel.: 0 825 02 30 30* (national)
Tel.: +33 (0) 1 45 14 85 00
(international)
info.fr@vwr.com
* 0,18 € TTC/min + prix appel

GERMANY

VWR International GmbH
Hilpertstraße 20a
D - 64295 Darmstadt
Tel.: 0800 702 00 07* (national)
Tel.: +49 (0) 6151 3972 0
(international)
info.de@vwr.com
*Freecall

HUNGARY

VWR International Kft.
Simon László u. 4.
4034 Debrecen
Tel.: +36 52 521130
info.hu@vwr.com

IRELAND

VWR International Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 (0) 1 88 22 222
sales.ie@vwr.com

ITALY

VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02 3320311
info.it@vwr.com

THE NETHERLANDS

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 (0) 20 4808 400
info.nl@vwr.com

NORWAY

VWR International AS
Brynsalleen 4,
0667 Oslo
Tel.: +47 22 90 00 00
info.no@vwr.com

POLAND

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 58 32 38 200
info.pl@vwr.com

PORTUGAL

VWR International –
Mat. de Laboratório, Soc.
Unipessoal, Lda
Edifício Ramazzotti
Avenida do Forte 6, P-1.09 e
P-1.10
2790-072 Carnaxide
Tel.: +351 21 3600 770
Info.pt@vwr.com

SPAIN

VWR International Eurolab S.L.U.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: +34 902 222 897
info.es@vwr.com

SWEDEN

VWR International AB
Fagerstagatan 18b
163 94 Stockholm
Tel.: +46 (0) 8 621 34 20
kundservice.se@vwr.com

SWITZERLAND

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: +41 (0) 44 745 13 13
info.ch@vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: +44 (0) 800 22 33 44
uksales@vwr.com

CHINA

VWR (Shanghai) Co., Ltd
Bld.No.1, No.3728 Jinke Rd,
Pudong New District
Shanghai, 201203- China
Tel.: 400 821 8006
info_china@vwr.com

INDIA

VWR Lab Products Private Limited
No.139, BDA Industrial Suburb,
6th Main, Tumkur Road, Peenya Post,
Bangalore, India – 560058
Tel.: +91 80 28078400
vwr_india@vwr.com

KOREA

Avantor Performance Materials Korea Ltd
2F ACE Gwanggyo Tower I, Daehak
4ro 17
Yeongtong-gu Suwon, Korea 16226
Tel.: +82 31 645 7256
saleskorea@avantorsciences.com

MIDDLE EAST & AFRICA

VWR International FZ-LLC
Office 203, DSP Lab Complex,
Dubai Science Park,
Dubai, United Arab Emirates
Tel.: +971 4 5573271
Info.mea@vwr.com

SINGAPORE

VWR Singapore Pte Ltd
The Metropolis
Tower 1, #05-03
9 North Buona Vista Drive
Singapore 138588
Tel.: +65 6505 0760
sales.sg@vwr.com

GO TO **VWR.COM** FOR THE LATEST NEWS, SPECIAL OFFERS AND DETAILS FROM YOUR LOCAL VWR SUPPORT TEAM