

Certificate of Conformity

VWR Universal pH Electrode, With 1 m Fixed Cable

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: 662-1763

Lot No: 0822

Type	Material	Connection	Electrolyte	ØxL (mm)	Pack
Flat tip, double junction (HDPE/Pellon)	Epoxy	BNC	Gel	12 x 150	1

pH range: 0 - 14

Temperature range: 10 to +100 °C

Response time: 95% \<5 seconds

Checked at 5 points during the manufacturing process and lot stamped to assure quality.

This document has been produced electronically and is valid without a signature

Date of issue: **June 2022**