

Specification

PBS tablets pH 7.4 (for 200 ml)

A9177

Product Code:	A9177
Product Name:	PBS tablets pH 7.4 (for 200 ml)
Specifications:	• Ready to use in minutes • pH (25°C): 7.4 ± 0.05 • Autoclaveable Composition: KCL: 2.7 mM NaCl: 140 mM Phosphate: 10 mM
WGK:	1
Storage:	RT
Origin:	Tablets for preparation of 200 ml phosphate buffered saline (PBS)
CS:	38221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340