



Certificate of Conformity

VWR Bulb Pipette, Two Mark, Class AS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Lot No.	Capacity ml	Division ml	Length mm	Tolerance	Pack
612-4775	GL1019	1	0.01	300	0.007 $\pm$ ml	3

Quality System Compliance: Products are manufactured under the ISO 9001:2008 certified plant & Calibrated in ISO 17025:2005 standard. Products are Inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release.

QC Testing: Representative production samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications:

The product is in conformance with the requirements of the following standards and regulations:

Technical standards for products: EN ISO 648

Calibration procedure: ISO 4787

Material:

The products are made from DIN ISO 3585:1998, TYPE-1, CLASS- A, Boro 3.3 Glass that meet ASTM E-438 DIN ISO 3585:1998, TYPE-1, CLASS- A Standard.

BPA statement: Not Applicable

Phthalate/DEHP Statement: Not Applicable

Heavy metals-RoHS: Not Applicable

Certificate of Conformity

REACH (SVHC) Statement: Not Applicable

BSE/TSE statement: The product is solely manufactured from Glass and is free from BSE and TSE. Glass is manufactured by using raw material of non-animal origin and does not contain meat product and does not come in contact with meat products.

Non-Pyrogenic: Not Applicable

DNase & RNase Free: Not Applicable

Cytotoxicity: Not Applicable

Latex Free statement: Not Applicable

Storage Conditions: Can be stored at room temperature in normal warehouse conditions

Shelf Life: There is no predefined shelf life assigned to the product.

Additional Information: The product is safe to use when used in proper way, as should be use with lab glassware. The purpose of the product is for use in laboratory only.

This document has been produced electronically and is valid without a signature

Date of issue: **May 2019**