

VWR International bv
Haasrode Researchpark 2020
Geldenaaksebaan 464
B-3001 Leuven
t 016 385 011

Certificate of Containment Testing

VWR Catalogue Number	521-2576
Description	Angle rotor RA 24/2 with aerosol tight lid
Country of Origin	Slovenia
Date of Issue	23 rd March 2022

Test summary:

Angle rotor RA 24/2 with aerosol tight lid was containment tested in a VWR centrifuge using Annex AA of standard EN 61010-2-020:2006.

Maximum speed of the rotor 18,000 rpm (= 31,150 rcf) was applied for 5 minutes.

The sealed rotor was shown to contain all contents.

This document has been produced electronically and is valid without a signature