

Molecular Biology

Mol Bio **Ethanol** **Nucleotide** **Ethanol** **RNase** **RNA** **DNA**
Molecular Biology
Nucleotide **RNase** **Ethanol** **DNA** **RNase** **Ethanol** **Mol Bio**
Nucleotide **Ethanol** **Mol Bio** **RNase** **DNA** **Nucleotide** **Ethanol** **Mol Bio** **RNase** **DNA**



460 Glennie Circle
King of Prussia, PA 19406
(800) 332-6647
www.deconlabs.com



Customer Service:
800-932-5000

MOLECULAR BIOLOGY

Molecular Biology Grade Ethanol

For Precipitation of Nucleic Acid

Molecular Biology Grade Ethanol effectively precipitates DNA and RNA from samples. 200 Proof (100% ethanol) can be used full strength or diluted with endonuclease-free water to desired dilution. Tested and shown to be absent of DNase and RNase, endonucleases that can degrade DNA and RNA. Available 24 x 1 Pint or sold as a single pint.



ELIMINase®

RNase/DNase Remover

ElimiNase® is a non-toxic, ready-to-use solution for instantly removing RNase, DNase, unwanted DNA, and radioactive contamination from labware and equipment. Safe and effective. Eliminates need for autoclaving or using DEPC.



PRODUCT GUIDE

PRODUCT	SIZE	VWR CAT. #	CASE SIZE
Molecular Biology Grade Ethanol <i>For Precipitation of Nucleic Acid</i>	1 Pint	71006-012	—
	1 Pint/24	71006-010	24 x 1 Pint
ELIMINase® <i>RNase/DNase Remover</i>	250 mL	10798-804	12 x 250 mL
	950 mL	10498-216	12 x 950 mL



460 Glennie Circle
King of Prussia, PA 19406
(800) 332-6647
www.deconlabs.com



Customer Service:
800-932-5000