

# Certificate of Conformity

|                                      |                                                                    |
|--------------------------------------|--------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>          | 732-3670                                                           |
| <b>Description</b>                   | <b>VWR® Gel-loading Pipette Round Tips</b><br>(Diameter – 0.57 mm) |
| <b>Lot Number</b>                    | 20492                                                              |
| <b>Date of Manufacture (yyyy-mm)</b> | 2020-12                                                            |
| <b>Country of Origin</b>             | Manufactured in USA                                                |
| <b>Date of Issue (yyyy-mm-dd)</b>    | 2021-08-12                                                         |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

## Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

|                               |                                                                           |
|-------------------------------|---------------------------------------------------------------------------|
| <b>Material</b>               | Polypropylene.<br>Volume of the tip – 200 µL<br>Length of the tip – 83 mm |
| <b>Sterilization</b>          | Not applicable (non-sterile).                                             |
| <b>ATP Assay</b>              | Not applicable.                                                           |
| <b>DNase &amp; RNase Free</b> | The above product is free of any detectable DNase/RNase contamination.    |

|                                |                                                                                                                                                                                                                                                                                                                                                                       |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).           |
| <b>Cytotoxicity</b>            | Testing has been conducted and have been shown to be non-toxic.                                                                                                                                                                                                                                                                                                       |
| <b>Non-Pyrogenic Statement</b> | The acceptance level for product is 0.06 EU/mL (LAL Gel Clot Method).                                                                                                                                                                                                                                                                                                 |
| <b>Latex Statement</b>         | The above product is not made from natural rubber latex.                                                                                                                                                                                                                                                                                                              |
| <b>BPA Statement</b>           | Bisphenol is not used in the manufacture of the raw material and are not expected to be present.                                                                                                                                                                                                                                                                      |
| <b>DEHP Statement</b>          | Please see RoHS information below.                                                                                                                                                                                                                                                                                                                                    |
| <b>RoHS</b>                    | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed. |
| <b>REACH Statement</b>         | No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.                                                                                                                                                      |
| <b>Storage Conditions</b>      | Store in room temperature.                                                                                                                                                                                                                                                                                                                                            |
| <b>Shelf Life</b>              | 5 years after manufacture.                                                                                                                                                                                                                                                                                                                                            |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.