



Product Specification Sheet

VWR PIPETTE TIPS, BULK

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Type	Capacity (µl)	Colour code	For	Packed	Pk
613-0259	UltraFine™, extended	0.1 - 10	-	-	Bulk	1.000
613-0276	UltraFine™, graduated, siliconised	0.1 - 10	-	-	Bulk	1.000
613-0364	UltraFine™, graduated	0.1 - 10	-	-	Bulk	1.000
613-0254	Bevel Point™, graduated	0.5 - 10	-	-	Bulk	1.000
613-0578	UltraFine™	0.5 - 20	-	Rainin® LTS™	Bulk	1.000
613-0636	Graduated	0.1 – 30	Transparent	VWR / Biohit / Gilson	Bulk	1.000
613-0252	Bevel Point™	1 - 100	-	Eppendorf / Brinkmann / Brand	Bulk	1.000
613-0239	Bevel Point™	1 - 200	Yellow	-	Bulk	1.000
613-0240	Bevel Point™	1 - 200	Yellow	-	Bulk	10.000
613-0246	Standard (yellow)	1 - 200	Yellow	-	Bulk	1.000
613-0277	Bevel Point™, siliconised	1 - 200	Green	-	Bulk	1.000
613-0287	Standard	1 - 200	Transparent	-	Bulk	1.000
613-0298	Bevel Point™, graduated	1 - 200	-	-	Bulk	1.000
613-0581	UltraFine™, graduated	1 - 200	-	Rainin® LTS™	Bulk	1.000
613-0890	Bevel Point™	1 - 200	Transparent	-	Bulk	1.000
613-0891	Bevel Point™, Low Binding	1 - 200	Transparent	-	Bulk	1.000
732-0500	UltraFine™, FlexTop™, extended	1 - 200	-	-	Bulk	1.000
732-0504	UltraFine™, FlexTop™, graduated	1 - 200	-	-	Bulk	1.000



Product Specification Sheet

VWR PIPETTE TIPS, BULK

Cat. No.	Type	Capacity (µl)	Colour code	For	Packed	Pk
613-0370	Low Binding, wide orifice	1 - 250	-	-	Bulk	1.000
613-0373	UltraFine™, Low Binding	1 - 250	-	-	Bulk	1.000
613-0892	UltraFine™	1 - 250	-	-	Bulk	1.000
736-0204	Bevel Point™, wide orifice	1 - 250	-	-	Bulk	1.000
613-0265	Bevel Point™	1 - 300	-	-	Bulk	1.000
613-0331	UltraFine™, FlexTop™, graduated	1 - 300	-	-	Bulk	1.000
613-0380	Bevel Point™, Low Binding	1 - 300	-	-	Bulk	1.000
525-0144	Bevel Point™, graduated	1 - 200	-	-	Bulk	10.000
613-0340	Bevel Point™	100 - 1000	-	-	Bulk	1.000
613-0343	Bevel Point™, graduated	100 - 1000	-	-	Bulk	1.000
613-0584	UltraFine™	100 - 1000	-	Rainin® LTS™	Bulk	1.000
613-0272	UltraFine™, FlexTop™, extended	100 - 1250	-	-	Bulk	1.000
613-0377	UltraFine™, Low Binding	100 - 1250	-	-	Bulk	1.000
613-0749	Wide orifice	100 - 1250	-	-	Bulk	1.000
613-0338	Macro, graduated (Gilson Pipetman)	1000 - 5000	Transparent	Gilson Pipetman	Bulk	250
613-0339	Macro, graduated (Biohit/Eppendorf)	1000 - 5000	Transparent	Biohit / Eppendorf	Bulk	250
613-0830	Macro	1000 - 5000	Transparent	VWR / Biohit / Gilson / Nichiryo / HTL / Labnet	Bulk	250
613-0831	Standard	5000 - 10000	-	Gilson Pipetman	Bulk	250

Material: 100% pure virgin medical grade PP



Product Specification Sheet

VWR PIPETTE TIPS, BULK

Other Characteristics:

- Made in the highest quality moulds, these tips can be used whenever the work demands a higher standard.
- Resins are pre-tested for metal contamination. Certificates are on file by lot number for all test results.
- A wider seal assures that tips will seal easily and consistently without leaking, on old as well as new pipettes
- Guaranteed free of detectable DNA, DNase/RNase, endotoxins and heavy metals
- Autoclavable

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **March 2013**

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium
<http://www.vwr.com>



Version: 2.0