

Certificate of Conformity

VWR Catalogue Number	441-0507
Description	Prime overhead stirrers, VWR®, 40 Ncm, 40 L Max. stirring capacity H ₂ O, 30000 mPas Max. Viscosity, 50 – 1000 min ⁻¹ speed range
Lot Number	QJ048976
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2022-07-11

Quality System Compliance

Products are manufactured under the **ISO 9001 & ISO 13485** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Instrument Body Material - Polypropylene and Impeller - SS316L.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Not Applicable.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.

DEHP Statement

Please see RoHS information below.

RoHS

No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.

REACH Statement

No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.

Storage Conditions

Permissible Ambient Temperature: 5-40° C.
Permissible Ambient Relative Humidity: 80%

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Annex

QC Report

Speed Calibration :

Sr. No.	Set RPM	Speed Accuracy	Actual RPM
1	50	45 - 55	53
2	500	495 - 505	497
3	1000	990 - 1010	995

Time Calibration :

Sr. No.	Set Value (Min)	Actual Time (Min)
1	01:00	00:01:00
2	05:00	00:05:00

Note: Deviation in time under ± 2 % is considered acceptable.

Equipment's Used for Calibration :

Sr. No.	Description	Certified By	Instrument Id No	Certificate No.	Calibration due Date
1	Digital Tachometer	PCC	Q-634367	20211123-3692-1	22-11-2022
2	Digital Stopwatch	PCC	SW-11	20211123-3692-4	22-11-2022

Note: 1) All instrument used for measurement checking are calibrated in NABL/NPL traceable lab.
2) Calibration traceable to national standards.
3) Result is related only to the item calibrated.