



Specification

1.00181.0002 Oxidase for the detection of cytochrome oxidase in microorganisms
Bactident®

Specification	
Cytochrome oxidase (<i>Pseudomonas aeruginosa</i> ATCC 10145 (WDCM 00024))	+
Cytochrome oxidase (<i>Pseudomonas aeruginosa</i> ATCC 27853 (WDCM 00025))	+
Cytochrome oxidase (<i>Pseudomonas aeruginosa</i> ATCC 9027 (WDCM 00026))	+
Cytochrome oxidase (<i>Pseudomonas fluorescens</i> ATCC 13525 (WDCM 00115))	+
Cytochrome oxidase (<i>Vibrio parahaemolyticus</i> NCTC 10885 (WDCM 00185))	+
Cytochrome oxidase (<i>Aeromonas hydrophila</i> ATCC 7966 (WDCM 00063))	+
Cytochrome oxidase (<i>Escherichia coli</i> ATCC 8739 (WDCM 00012))	-
Cytochrome oxidase (<i>Escherichia coli</i> ATCC 25922 (WDCM 00013))	-
Cytochrome oxidase (<i>Escherichia coli</i> ATCC 11775 (WDCM 00090))	-
Cytochrome oxidase (<i>Escherichia coli</i> DSM 19683 (WDCM 00179))	-
Cytochrome oxidase (<i>Cronobacter sakazakii</i> ATCC 29544 (WDCM 00214))	-
Cytochrome oxidase (<i>Cronobacter muytjensis</i> ATCC 51329 (WDCM 00213))	-
Cytochrome oxidase (<i>Serratia marcescens</i> ATCC 14756)	-

Stefanie Fischer
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.