

Product Specification

Test	Specification
Assay	49.5 to 50.5 %
Ammonium Hydroxide Precipitate	<= 0.01 %
Chloride (Cl)	<= 0.01 %
Heavy Metals (as Ag)	<= 0.001 %
Potassium (K)	<= 0.01 %
Sodium Carbonate (Na ₂ CO ₃)	<= 0.5 %
Sulphate (SO ₄)	<= 0.002 %
Trace Impurities (in ppm)	
Iron (Fe)	<= 5
Nickel (Ni)	<= 5
Nitrogen Compounds (as N)	<= 5
Phosphate (PO ₄)	<= 5

For Laboratory, Research or Manufacturing use.

Storage conditions: IMPORTANT: Material will freeze if stored below 18°C.



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000