

1121 5th Street NW
Washington, DC 20001
www.edvotek.com

Safety Data Sheet

Version 4
Revision Date 05/14/2025

1. PRODUCT AND COMPANY IDENTIFICATION

Product name: Colored Dye (Red,
Blue, Purple, Yellow, Navy)
Product Number:
Brand: Edvotek
Supplier: Edvotek
1121 5th Street NW
Washington, DC
20001

Telephone: 1-800-338-6835
or 202-370-1500
Fax: 202-370-1501
Emergency Phone #: 202-370-1500

2. HAZARDS IDENTIFICATION

Emergency Overview

OSHA Hazards:

No known OSHA hazards.

Not a dangerous substance or mixture according to the Globally Harmonized System (GHS).

HMIS Classification:

Health hazard: 0
Flammability: 0
Physical hazard: 0

NFPA Rating:

Health hazard: 0
Fire: 0
Reactivity Hazard: 0

Potential Health Effects:

Inhalation: May be harmful if inhaled.
Skin: May be harmful if absorbed through skin. May cause skin irritation.
Eyes: May cause eye irritation.
Ingestion: May be harmful if swallowed.

3. COMPOSITION/INFORMATION ON INGREDIENTS

No ingredients are hazardous according to OSHA criteria.

4. FIRST AID MEASURES

General information

No special measures required.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult doctor in case of complaints.

In case of skin contact

Immediately wash off with soap and plenty of water. Get medical attention if symptoms develop.

In case of eye contact

Flush eyes with water as a precaution. Get medical attention immediately if irritation persists.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Get medical attention if symptoms develop.

5. FIREFIGHTING MEASURES

Conditions of flammability

Not flammable or combustible.

Suitable extinguishing media

Use water spray, alcohol resistant foam, dry chemical or carbon dioxide.

Special protective equipment for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

Hazardous combustion products

No known products

6. ACCIDENTAL RELEASE MEASURES

Personal precautions

Prevent skin/eye contact. Use personal protective equipment as needed.

Environmental precautions

Dilute with plenty of water. Do not allow to enter sewers/ surface or ground water.

Methods and materials for containment and cleaning up

Clean up spills immediately. Use absorbent materials if needed and dispose into a chemical waste container.

7. HANDLING AND STORAGE

Precautions for safe handling

No special measures necessary. Good laboratory technique should be used when handling.

Conditions for safe storage

No special measures necessary. Store at room temperature.

Incompatible Products

No known products

8. EXPOSURE CONTROLS/ PERSONAL PROTECTION

Personal protective equipment

Respiratory protection

Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

Wear protective gloves. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices.

Eye protection

Avoid eye contact with this material. Wear chemical safety goggles.

Skin and body protection

Use chemical gloves as needed. Use clean protective body covering clothing as needed to minimize contact with clothing and skin.

Hygiene measures

Handle in accordance with good industrial hygiene and safety practice. Wash and dry hands at the end of workday.

9. PHYSICAL AND CHEMICAL PROPERTIES**Appearance**

Form	liquid
Color	red / blue/ yellow
Odor	no data available

Safety data

pH	no data available
Melting point/freezing point	no data available
Boiling point	no data available
Flash point	no data available
Ignition temperature	no data available
Autoignition Temperature	no data available
Lower explosion limit	no data available
Upper explosion limit	no data available
Vapour pressure	no data available
Density (water=1)	no data available
Water solubility	no data available
Partition coefficient:	
n-octanol/water	no data available
Relative vapour Density	no data available
Odour Threshold	no data available
Evaporation rate	no data available

10. STABILITY AND REACTIVITY

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

No further relevant data available.

Conditions to avoid

No further relevant data available.

Materials to avoid

No further relevant data available.

Hazardous decomposition products

No dangerous decomposition products known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity LD/LC50 values

no data available

Serious eye damage/eye irritation

no data available

Respiratory or skin sensitization

no data available

Germ cell mutagenicity

no data available

Skin corrosion/irritation

no data available

Carcinogenicity

No component of this product are listed by IARC, ACGIH, NTP and OSHA.

Reproductive toxicity

no data available

Teratogenicity

no data available

Aspiration hazard

no data available

Synergistic effects

no data available

12. ECOLOGICAL INFORMATION

Toxicity

no data available

Persistence and degradability

no data available

Bioaccumulative potential

no data available

Mobility in soil

no data available

PBT and vPvB assessment

no data available

Other adverse effects

no data available

13. DISPOSAL CONSIDERATIONS

Product

Dispose of unused product, residues, and containers according to local, regional, state, and national regulations.

Contaminated packaging

Dispose of as unused product.

14. TRANSPORT INFORMATION

DOT (US)

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

15. REGULATORY INFORMATION

OSHA Hazards

No known OSHA hazards.

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards.

16. OTHER INFORMATION

For educational use only. Not intended for human or animal diagnostic or therapeutic uses.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.



1121 5th Street NW
Washington, DC 20001
www.edvotek.com

SAFETY DATA SHEET

Version 1
Revision Date 05/14/2025

1. PRODUCT AND COMPANY IDENTIFICATION

Product name:	Clear Blue Fluorescent Dye Fluorescent Agent	Telephone:	1-800-338-6835 or 202-370-1500
Product Number:		Fax:	202-370-1501
Brand:	Edvotek	Emergency Phone #:	202-370-1500
Supplier:	Edvotek 1121 5 th Street NW Washington, DC 20001	(For both supplier and Manufacturer)	

2. HAZARDS IDENTIFICATION

Emergency Overview OSHA Hazards:

Not Established Triethanolamine ACGIH:5.00 mg/m³
TWA 10-15% Urea OSHA: Not Established ACGIH: Not
Established OSHA 10-15%

NFPA Rating:

Health hazard: 0
Fire: 0
Reactivity Hazard: 0

HMIS Classification:

Health hazard: 0
Flammability: 0
Physical hazard: 0

Potential Health Effects:

Skin: May cause skin irritation.
Eyes: May cause eye irritation.
Ingestion: May be harmful if swallowed.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Emergency Overview: Color: yellow; Form: Liquid;
Odor: Ammonia; May cause eye, skin, and respiratory
tract irritation.

injury to the mouth throat and digestive tract with
abdominal pain, burning, vomiting, and diarrhea.

Acute Eye Contact: Due to the low vapor pressure
inhalation of this material would be unlikely at room
temperature. If heated vapors or mists may cause
respiratory tract irritation with symptoms of coughing,
runny nose and sore throat.

Chronic Effects of Exposure: No applicable
information was found concerning any adverse chronic
health effects.

Acute Skin Contact: May cause irritation to the skin
with symptoms of drying, itching and reddening.

Carcinogenicity: The components of this product are
not listed by NTP, IARC, or regulated as a carcinogen by
OSHA.

Acute Ingestion: May cause nausea and produce local

Medical Conditions Aggravated by Exposure: Person
with preexisting eye, skin, or respiratory tract disorders
may be more susceptible to the effects of this product.

4. FIRST AID MEASURES

Eye Contact: Promptly flush eyes with large amounts of water for at least 15 minutes. Get prompt medical attention.

Skin Contact: If large amounts then immediate flush with large amounts of water use soap if available. Remove contaminated clothing including shoes after flushing has begun. Call physician if irritation persists.

Inhalation: If inhaled immediately remove the affected victim from exposure. Administer artificial respiration if breathing is stopped. Keep at rest. Call for medical attention.

Ingestion: If swallowed call for medical attention

5. FIREFIGHTING MEASURES

FlashPoint: >200 F

Extinguishing Media: Water; dry chemical; carbon dioxide; foam

Special Fire Fighting Procedures: Firefighters should be equipped with self-contained breathing apparatus to protect against potentially toxic and irritation fumes.

6. ACCIDENTAL RELEASE MEASURES

Spill or Leak: Utilize recommended protective clothing. Spills should be taken up with an appropriate absorbent and placed in approved containers. Spill area can be washed with water. Collect wash water for approved disposal.

7. HANDLING AND STORAGE

Shelf Life: 12 Months.

Storage Temperature: 40 F Minimum

Storage and Handling: Stir before using. Avoid breathing mist vapors. Avoid contact with eyes and skin (in concentrated form). Wash thoroughly after handling. Store in a dry place away from excessive heat, in original or similar waterproof containers. Re seal containers immediately after use. Keep from freezing.

8. EXPOSURE CONTROLS/ PERSONAL PROTECTION

Eye Protection Requirements: Chemical Safety goggles.

conditions.

Personal Protection: Rubber gloves and long sleeved shirts and pants to minimize skin contact. Wash hands before eating and drinking. NO special ventilation or respirator requirements when used under normal

Additional Protective Measures: Emergency showers and eye wash stations should be available. Educate and train employees in the safe use and handling of this product.

9. PHYSICAL AND CHEMICAL PROPERTIES

Specific Gravity AF: 1.1

Solubility in Water: Miscible

Color: Yellow

Physical Form: Liquid

VOC by Weight: Approx. 2% EPA method 24/24a

pH: 10-11

Freezing/melting point AF: 23/212 F

10. STABILITY AND REACTIVITY

Stability: Stable

Oxidizing and reducing agents may destroy color.

Conditions to Avoid instability: Not applicable

Hazardous Decomposition Products: CO, CO₂ oxides of nitrogen oxides of sulfur and other potentially toxic fumes.

Hazardous Polymerization: Will not occur

Materials and Conditions to Avoid incompatibility:

11. TOXICOLOGICAL INFORMATION

Oral LD50: 14,530 mg/kg (rat)

Dermal LD50: >2000 mg/kg (rabbit)

12. ECOLOGICAL INFORMATION

Toxicity
no data available

PBT and vPvB assessment
no data available

Persistence and degradability
no data available

Other adverse effects
no data available

Bioaccumulative potential
no data available

Mobility in soil
no data available

13. DISPOSAL CONSIDERATIONS

Product

Dispose of unused product, residues, and containers according to local, regional, state, and national regulations.

Contaminated packaging

Dispose of as unused product.

14. TRANSPORT INFORMATION

DOT (US)	Not dangerous goods
IMDG	Not dangerous goods
IATA	Not dangerous goods

15. REGULATORY INFORMATION

OSHA Hazards

No known OSHA hazards

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards.

16. OTHER INFORMATION

For educational use only. Not intended for human or animal diagnostic or therapeutic uses.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein.

THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.



1121 5th Street NW
Washington, DC 20001
www.edvotek.com

SAFETY DATA SHEET

Version 1
Revision Date 05/14/2025

1. PRODUCT AND COMPANY IDENTIFICATION

Product name:	Fluorescent Yellow Green UV Tracer	Telephone:	1-800-338-6835 or 202-370-1500
Product Number:		Fax:	202-370-1501
Brand:	Edvotek	Emergency Phone #:	202-370-1500
Supplier:	Edvotek 1121 5 th Street NW Washington, DC 20001		(For both supplier and Manufacturer)

2. HAZARDS IDENTIFICATION

Emergency Overview OSHA Hazards:

Not Established Triethanolamine ACGIH: 5.00 mg/m³
TWA 10-15% Urea OSHA: Not Established ACGIH: Not
Established OSHA 10-15%

NFPA Rating:

Health hazard: 0
Fire: 0
Reactivity Hazard: 0

HMIS Classification:

Health hazard: 0
Flammability: 0
Physical hazard: 0

Potential Health Effects:

Skin: May cause skin irritation.
Eyes: May cause eye irritation.
Ingestion: May be harmful if swallowed.

3. COMPOSITION/INFORMATION ON INGREDIENTS

This product is not hazardous according to OSHA 29 CFT 1910.1200. Components not listed are not hazardous or are below reportable limits.

4. FIRST AID MEASURES

Eye Contact: Promptly flush eyes with plenty of water for at least 15 minutes. If irritation persists, get prompt medical attention.

Inhalation: Remove to fresh air. If breathing difficult, administer oxygen and seek medical attention.

Skin Contact: Wash thoroughly with water and soap. If irritation persists, get prompt medical attention.

Ingestion: Rinse mouth. DO NOT induce vomiting. Drink plenty of water. Seek medical attention if large quantities were ingested or if nausea occurs.

5. FIREFIGHTING MEASURES

Suitable Extinguishing Media: Water spray (fog); dry chemical; carbon dioxide (CO₂)

Special Fire Fighting Procedures: Product is not flammable. Burning/combustion may produce oxides of carbon and nitrogen.

Protective Equipment and Precautions for

Firefighters: Wear Self-contained breathing apparatus, MSHA/NIOSH (approved or equivalent) and full protective gear.

6. ACCIDENTAL RELEASE MEASURES

Spill or Leak: Use personal protective equipment as recommended in section 8. Prevent from entering into the soil, ditches, sewers, waterways, groundwater. Contain spill to prevent further leakage. Sweep up and collect into a suitable container for disposal. Flush area with water.

7. HANDLING AND STORAGE

Safe handling: Handle in accordance with good industrial hygiene and safety practices. Use personal protection recommended in section 8. Avoid contact with skin, eyes, or clothing. Avoid breathing dusts. Contaminated clothing should not be allowed out of the workplace.

Storage conditions: Keep container tightly closed and store in a cool, dry and well-ventilated area. Keep from freezing.

Incompatible Materials: Acids

8. EXPOSURE CONTROLS/ PERSONAL PROTECTION

This product does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

Engineering Controls: Ensure adequate ventilation, especially in confined areas. Eyewash stations. Showers.

Eye Protection Requirements: Chemical Safety goggles.

Personal Protection: Rubber gloves and suitable protective clothing. NO special ventilation or respirator requirements when used under normal conditions.

Hygiene Measures: Handle in accordance with good industrial hygiene and safety practices..

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance

Form	liquid
Colour	yellow/green liquid
Odor	No odor

Safety data

pH	>8.0	Vapour pressure	no data available
Melting point/freezing point	~ 32°F	Vapour Density (water=1)	0.6
Boiling point	~ 212°F	Water solubility	Highly soluble
Flash point	no data available	Partition coefficient:	
Ignition temperature	no data available	n-octanol/water	no data available
Autoignition Temperature	no data available	Relative vapour Density	no data available
Lower explosion limit	no data available	Odour Threshold	no data available
Upper explosion limit	no data available	Evaporation rate	1.8

10. STABILITY AND REACTIVITY

Stability: Stable

Materials and Conditions to Avoid incompatibility:
Acids

Conditions to Avoid instability: Not applicable

Hazardous Decomposition Products: Oxides of carbon and nitrogen.

Hazardous Polymerization: Will not occur

11. TOXICOLOGICAL INFORMATION

Inhalation: Avoid breathing vapors or mists.

Ingestion: Do not ingest.

Skin contact: may cause an allergic skin reaction.

Eye contact: avoid contact with eyes.

Acute toxicity LD/LC50 values

no data available

Serious eye damage/eye irritation

no data available

Respiratory or skin sensitization

no data available

Germ cell mutagenicity

no data available

Skin corrosion/irritation

no data available

Carcinogenicity

No component of this product are listed by IARC, ACGIH, NTP and OSHA

Reproductive toxicity

no data available

Teratogenicity

no data available

Aspiration hazard

no data available

Synergistic effects

no data available

12. ECOLOGICAL INFORMATION

Toxicity

no data available

PBT and vPvB assessment

no data available

Persistence and degradability

no data available

Other adverse effects

no data available

Bioaccumulative potential

no data available

Mobility in soil

no data available

13. DISPOSAL CONSIDERATIONS

Product

Dispose of unused product, residues, and containers according to local, regional, state, and national regulations.

Contaminated packaging

Dispose of as unused product.

14. TRANSPORT INFORMATION

DOT (US)

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

15. REGULATORY INFORMATION

OSHA Hazards

No known OSHA hazards

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the reporting levels established by SARA Title III, Section 313.

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 311/312 Hazards


No SARA Hazards.

16. OTHER INFORMATION

For educational use only. Not intended for human or animal diagnostic or therapeutic uses.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein.

THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.



1121 5th Street NW
Washington, DC 20001
www.edvotek.com

Safety Data Sheet

Version 3
Revision Date 05/14/2025

1. PRODUCT AND COMPANY IDENTIFICATION

Product name:	Simulated Blood Sample	Telephone:	1-800-338-6835 or 202-370-1500
Product Number:		Fax:	202-370-1501
Brand:	Edvotek	Emergency Phone #:	202-370-1500
Supplier:	Edvotek 1121 5 th Street NW Washington, DC 20001		

2. HAZARDS IDENTIFICATION

Emergency Overview

OSHA Hazards:

No known OSHA hazards.

Not a dangerous substance or mixture according to the Globally Harmonized System (GHS).

HMIS Classification:

Health hazard: 0
Flammability: 0
Physical hazard: 0

NFPA Rating:

Health hazard: 0
Fire: 0
Reactivity Hazard: 0

Potential Health Effects:

Inhalation: May be harmful if inhaled.
Skin: May be harmful if absorbed through skin. May cause skin irritation.
Eyes: May cause eye irritation.
Ingestion: May be harmful if swallowed.

3. COMPOSITION/INFORMATION ON INGREDIENTS

No ingredients are hazardous according to OSHA criteria.

4. FIRST AID MEASURES

General information

No special measures required.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult doctor in case of complaints.

In case of skin contact

Immediately wash off with soap and plenty of water. Get medical attention if symptoms develop.

In case of eye contact

Flush eyes with water as a precaution. Get medical attention immediately if irritation persists.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Get medical attention if symptoms develop.

5. FIREFIGHTING MEASURES

Conditions of flammability

Not flammable or combustible.

Suitable extinguishing media

Use water spray, alcohol resistant foam, dry chemical or carbon dioxide.

Special protective equipment for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

Hazardous combustion products

No known products

6. ACCIDENTAL RELEASE MEASURES

Personal precautions

Prevent skin/eye contact. Use personal protective equipment as needed.

Environmental precautions

Dilute with plenty of water. Do not allow to enter sewers/ surface or ground water.

Methods and materials for containment and cleaning up

Clean up spills immediately. Use absorbent materials if needed and dispose into a chemical waste container.

7. HANDLING AND STORAGE

Precautions for safe handling

No special measures necessary. Good laboratory technique should be used when handling.

Conditions for safe storage

No special measures necessary. Store at room temperature.

Incompatible Products

No known products

8. EXPOSURE CONTROLS/ PERSONAL PROTECTION

Personal protective equipment

Respiratory protection

Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

Wear protective gloves. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices.

Eye protection

Avoid eye contact with this material. Wear chemical safety goggles.

Skin and body protection

Use chemical gloves as needed. Use clean protective body covering clothing as needed to minimize contact with clothing and skin.

Hygiene measures

Handle in accordance with good industrial hygiene and safety practice. Wash and dry hands at the end of workday.

9. PHYSICAL AND CHEMICAL PROPERTIES**Appearance**

Form liquid
Color Red
Odor no data available

Safety data

pH	no data available
Melting point/freezing point	no data available
Boiling point	no data available
Flash point	no data available
Ignition temperature	no data available
Autoignition Temperature	no data available
Lower explosion limit	no data available
Upper explosion limit	no data available
Vapor pressure	no data available
Density (water=1)	no data available
Water solubility	no data available
Partition coefficient:	
n-octanol/water	no data available
Relative vapour Density	no data available
Odour Threshold	no data available
Evaporation rate	no data available

10. STABILITY AND REACTIVITY

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

No further relevant data available.

Conditions to avoid

No further relevant data available.

Materials to avoid

No further relevant data available.

Hazardous decomposition products

No dangerous decomposition products known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity LD/LC50 values

no data available

Serious eye damage/eye irritation

no data available

Respiratory or skin sensitization

no data available

Germ cell mutagenicity

no data available

Skin corrosion/irritation

no data available

Carcinogenicity

No components of this product are listed by IARC, ACGIH, NTP and OSHA.

Reproductive toxicity

no data available

Teratogenicity

no data available

Aspiration hazard

no data available

Synergistic effects

no data available

12. ECOLOGICAL INFORMATION**Toxicity**

no data available

Persistence and degradability

no data available

Bioaccumulative potential

no data available

Mobility in soil

no data available

PBT and vPvB assessment

no data available

Other adverse effects

no data available

13. DISPOSAL CONSIDERATIONS**Product**

Dispose of unused product, residues, and containers according to local, regional, state, and national regulations.

Contaminated packaging

Dispose of as unused product.

14. TRANSPORT INFORMATION**DOT (US)**

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

15. REGULATORY INFORMATION**OSHA Hazards**

No known OSHA hazards.

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards.

16. OTHER INFORMATION

For educational use only. Not intended for human or animal diagnostic or therapeutic uses.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.