



MP Biomedicals, LLC

29525 Fountain Parkway

Solon OH 44139

United States

800.854.0530

www.mpbio.com

Certificate of Analysis

Description: Minimum Essential Medium Eagle With Earle's Salts and 2.0 g/L Sodium Bicarbonate Without L-Glutamine, L-Cysteine, L-Cystine, and L-Methionine

Catalog number: 091641454

QI Request No: 00007123

Related Catalog number:

Lot number: C2561

Quality Level: N/A

Recommended Retest Date (MM/DD/YYYY): 02-28-2025

Formula: N/A

Formula Weight: N/A

CAS #: N/A

Storage: 2 - 8°C

Storage Conditions: N/A

Additional Product information

Additional Lot Information

Test Description	Specification	Result	Comments
Identity	PASS	PASS	
Endotoxin (EU/ml)	<=0.25	<0.10	
Functionality Test (Passes)	PASS	PASS	
Osmolality (mM/kg)	260 - 320	282	
pH	7.0 - 7.4	7.4	
Physical Description (Clear Red Solution)	PASS	PASS	Clear Red Solution
Sterility (Passes)	PASS	PASS	

05-19-2023 Priti Faldu

MP Biomedicals, LLC.

This is an electronically generated document

Mail to: biotech@mpbio.com