



Product Specification Sheet

VWR SQUARE GLASS CUVETTES, LIGHT PATH 20 mm

VWR Cat.No.	VWRI634-6014
Product Name	SQUARE.GLASS CUVETTE (20mm)
Serial Number	634-6014
Material	B270 Optical Glass
Material Thickness	1.25±0.1 (mm)
Size	22.5*12.5*45 L*W*H (mm)
Optical Path Length	20±0.05 (mm)
Optical Path Width	10±0.1 (mm)
Cubic Measure	7 (ml)
Transmittance	≥79% (at 700nm)
Package	2pcs/box
Deviations in Matching	≤0.5%

VWR International bvba/sprl, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>

Date of issue: Dec.03rd 2012