


VWR®

Universal Low Retention Pipet Tips

01. CLEAN DELIVERY

02. PERFECT FIT





Compatible with universal pipettors including Gilson®, Eppendorf®, Biohit®, and more!

Secure Seal

1 FLEXIBLE PROXIMAL END

- Gives a consistently secure seal with the widest range of industry-leading pipettors
- Requires less insertion force than competitors' to make a secure seal, reducing the risk of repetitive stress injury
- Tips don't fall off your multichannel pipettor
- Makes every pipettor more ergonomic

Clean Delivery

2 LOW RETENTION SURFACE

- Maximizes sample recovery
- Improves accuracy across your experiments

3 MINIMIZES SURFACE TENSION

- Eliminates hanging droplet formation
- Saves time by reducing the need for tip touch off
- Increases accuracy and precision

VWR® Universal Pipet Tips, Low Retention, Filtered (960 Tips/Pack)

Volume	2 µl	10 µl	10 µl, XL	20 µl	100 µl	200 µl	300 µl	1000 µl	1250 µl
Sterilized	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Racked	732-3709	732-3630	732-3631	732-3632	732-3633	732-3634	732-3636	732-3637	732-3638

VWR® Universal Pipet Tips, Low Retention, Non-Filtered (960 Tips/Pack)

Volume	10 µl		10 µl, XL		200 µl		250 µl		300 µl		1000 µl		1250 µl	
Sterilized	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes
Racked	613-6453	613-6462	613-6454	613-6463	613-6456	613-6456	613-6458	613-6468	613-6459	613-6469	613-6460	613-6470	613-6461	613-6471
Reload	613-6433	613-6440	613-6434	613-6441	613-6435	613-6442	613-6436	–	613-6437	613-6443	613-6438	–	613-6439	613-6444

Empty Racks

Volume	10 - 300 µl	1000 - 1250 µl
Cat. No.	613-6445	613-6446