



VWR® 2000mL Erlenmeyer Flask

Description

2000 mL, PC, Non-treated, Plug Seal Cap or Vent Cap, Sterile

Purpose

Ideal for shaker culture applications and storage


Materials

Flask: PC (Polycarbonate)

Cap: HDPE (High-density Polyethylene) Color: Green

Membrane: 0.22 µm Hydrophobic PTFE (only use for Vent Cap)

Dimensions (Unit: mm)



Features

- Temperature range for PC Flasks: -80°C to +121°C
- PC provides excellent optical clarity and mechanical strength
- Unlike glass, flasks will not break when dropped
- Molded-in graduations for accuracy
- Available with plug seal cap or vent cap
- Vent cap option for continuous gas exchange while ensuring sterility and preventing leakage
- Triple bags, individually packaged for VWR 2L-3L-5L PC Erlenmeyer Flasks
- Non-cytotoxic
- DNase/RNase free and non-pyrogenic
- Sterilized by irradiation SAL 10⁻⁶ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system



STERILE

R

VWR® 2000mL Erlenmeyer Flask

VWR EU Cat. No.	Capacity	Flask Material	Cap Style	Cap Color	Sterilization	Packaging Configuration
214-0476	2000 mL	PC	Plug seal cap	Green	Yes	1/PK, 6 PK/CS
214-0477	2000 mL	PC	Vent cap	Green	Yes	1/PK, 6 PK/CS
214-0488	Cap		Plug seal cap	Green	Yes	1/PK, 24 PK/CS
214-0489	Cap		Vent cap	Green	Yes	1/PK, 24 PK/CS



VWR® 3000mL Erlenmeyer Flask

Description

3000 mL, PC, Non-treated, Plug Seal Cap or Vent Cap, Sterile

Purpose

Ideal for shaker culture applications and storage

Materials


Flask: PC (Polycarbonate)

Cap: HDPE (High-density Polyethylene)

Color: Green

Membrane: 0.22 µm Hydrophobic PTFE (only use for Vent Cap)

Dimensions (Unit: mm)



Features

- Temperature range for PC Flasks: -80°C to +121°C
- PC provides excellent optical clarity and mechanical strength
- Unlike glass, flasks will not break when dropped
- Molded-in graduations for accuracy
- Available with plug seal cap or vent cap
- Vent cap option for continuous gas exchange while ensuring sterility and preventing leakage
- Triple bags, individually packaged for VWR 2L-3L-5L PC Erlenmeyer Flasks
- Non-cytotoxic
- DNase/RNase free and non-pyrogenic
- Sterilized by irradiation SAL 10^{-6} (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system



STERILE

R

VWR® 3000mL Erlenmeyer Flask

VWR EU Cat. No.	Capacity	Flask Material	Cap Style	Cap Color	Sterilization	Packaging Configuration
214-0478	3000 mL	PC	Plug seal cap	Green	Yes	1/PK, 4 PK/CS
214-0479	3000 mL	PC	Vent cap	Green	Yes	1/PK, 4 PK/CS
214-0490	Cap		Plug seal cap	Green	Yes	1/PK, 24 PK/CS
214-0491	Cap		Vent cap	Green	Yes	1/PK, 24 PK/CS



VWR® 5000mL Erlenmeyer Flask

Description

5000 mL, PC, Non-treated, Plug Seal Cap or Vent Cap, Sterile

Purpose

Ideal for shaker culture applications and storage

Materials


Flask: PC (Polycarbonate)

Cap: HDPE (High-density Polyethylene)

Color: Green

Membrane: 0.22 µm Hydrophobic PTFE (only use for Vent Cap)

Dimensions (Unit: mm)



Features

- Temperature range for PC Flasks: -80°C to +121°C
- PC provides excellent optical clarity and mechanical strength
- Unlike glass, flasks will not break when dropped
- Molded-in graduations for accuracy
- Available with plug seal cap or vent cap
- Vent cap option for continuous gas exchange while ensuring sterility and preventing leakage
- Triple bags, individually packaged for VWR 2L-3L-5L PC Erlenmeyer Flasks
- Non-cytotoxic
- DNase/RNase free and non-pyrogenic
- Sterilized by irradiation SAL 10^{-6} (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system



STERILE

R

VWR® 5000mL Erlenmeyer Flask

VWR EU Cat. No.	Capacity	Flask Material	Cap Style	Cap Color	Sterilization	Packaging Configuration
214-0480	5000 mL	PC	Plug seal cap	Green	Yes	1/PK, 4 PK/CS
214-0481	5000 mL	PC	Vent cap	Green	Yes	1/PK, 4 PK/CS
214-0492	Cap		Plug seal cap	Green	Yes	1/PK, 24 PK/CS
214-0493	Cap		Vent cap	Green	Yes	1/PK, 24 PK/CS