

Sterile Millex® syringe filters you can trust.

Research Applications:

Sterile filtration of tissue culture media and additives / Buffers / DMSO / Biological solutions

Membrane	Pore Size (µm)	Diameter (mm)	Process Volume (hold-up)	Housing, sterilization method	Connections (inlet outlet)	10 PK	25 PK	50 PK	100 PK	250 PK
PES Millipore Express® Plus Membrane Fast flow and low binding for cell culture media preparation	0.22	25	100 mL (< 100 µL)	PVC, EO	FLL MLS			SLMPR25SS		
		33	200 mL (< 100 µL)	Modified acrylic, RS	FLL MLS			SLGPR33RS		SLGPR33RB
		50	4-8 L	Acrylic, RS	HB HB	SLGP05010				
					HB HB with filling bell	SLGPB5010				
	0.45	33	200 mL (< 100 µL)	Modified acrylic, RS	FLL MLS			SLHPR33RS		SLHPR33RB
PVDF Durapore® Membrane Lowest binding for protein rich solutions	0.1	33	100 mL (< 100 µL)	Modified acrylic, RS	FLL MLS			SLVVR33RS		
	0.22	4	1 mL (< 10 µL)	HDPE, EO	FLL Male stepped				SLGV004SL	
		13	10 mL (< 25 µL)	HDPE, EO	FLL MLS				SLGVR13SL	
		25	100 mL (< 100 µL)	PVC, EO	Vented FLL MLS			SLGVVR255F		
		33	100 mL (< 100 µL)	Modified acrylic, RS	FLL MLS			SLGVR33RS		SLGVR33RB
	0.45	4	1 mL (< 10 µL)	HDPE, EO	FLL Male stepped				SLHV004SL	
		13	10 mL (< 25 µL)	HDPE, EO	FLL MLS				SLHVR13SL	
		33	100 mL (< 100 µL)	Modified acrylic, RS	FLL MLS			SLHVR33RS		SLHVR33RB
	5.0	25	100 mL (< 100 µL)	PVC, EO	FLL MLS			SLSV0R25LS		
PVDF Durapel™ Membrane Superhydrophobic for air venting and gas filtration	0.22	25	N/A	PVC, EO	FLL MLL			SLGVS25US		
MCE MF-Millipore™ Membrane Most referenced general purpose membrane	0.22	25	100 mL (< 100 µL)	PVC, EO	Vented FLL MLS			SLGSVR255F		
		33	100 mL (< 100 µL)	Modified acrylic, EO	FLL MLS			SLGSR33SS		SLGSR33SB
	0.45	33	100 mL (< 100 µL)	Modified acrylic, EO	FLL MLS			SLHAR33SS		SLHAR33SB
	0.8	33	100 mL (< 100 µL)	Modified acrylic, EO	FLL MLS			SLAAR33SS		SLAAR33SB
Hydrophilic PTFE Membrane Broad chemical compatibility	0.2	13	10 mL (< 25 µL)	HDPE, EO	FLL MLS				SLLG013SL	
		25	100 mL (< 100 µL)	HDPE, EO	FLL MLS			SLLG025SS		
PTFE Fluoropore™ Membrane Hydrophobic for gas and alcohol filtration	0.2	25	100 mL (< 100 µL)	PVC, EO	FLL MLS			SLFGR25LS		
					FLL MLL			SLFGL25BS		
					FLL MLS with Needle		SLFGNR25VS			

Medical Applications*: CE marked and/or FDA registered for medical intended use.

Membrane	Pore Size (µm)	Diameter (mm)	Process Volume (hold-up)	Housing, sterilization method	Connections (inlet outlet)	CE	FDA (US)	25 PK	50 PK	100 PK	250 PK
PES Millipore Express® Plus Membrane Fast flow and low binding	0.22	25	100 mL (< 100 µL)	PVC, EO	FLL MLS	✓	✓		SLMP025SS		
					FLL MLL	✓	✓		SLMPL25SS		
		33	200 mL (< 100 µL)	Modified acrylic, RS	FLL MLS	✓	✓		SLGP033RS		SLGP033RB
	0.45	33	200 mL (< 100 µL)	Modified acrylic, RS	FLL MLS	✓	✓		SLGPM33RS		
						✓	✓		SLHP033RS		SLHP033RB
PVDF Durapore® Membrane Lowest protein binding	0.1	33	100 mL (< 100 µL)	Modified acrylic, RS	FLL MLS	✓	✓		SLV033RS		
						✓	✓		SLVVM33RS		
	0.22	13	10 mL (< 25 µL)	HDPE, EO	FLL MLS	✓	✓			SLGV013SL	
		25	100 mL (< 100 µL)	PVC, EO	Vented FLL MLS	✓	✓		SLGVV255F		
					Vented FLL MLL	✓	✓		SLGV0250S		
					FLL MLS	✓	✓		SLGV033RS		SLGV033RB
		33	100 mL (< 100 µL)	Modified acrylic, RS	FLL MLS	✓	✓		SLGVM33RS		
	0.45	13	10 mL (< 25 µL)	HDPE, EO	FLL MLS	✓	✓			SLHV013SL	
		33	100 mL (< 100 µL)	Modified acrylic, RS	FLL MLS	✓	✓		SLHV033RS		SLHV033RB
						✓	✓		SLHVM33RS		
	5.0	25	100 mL (< 100 µL)	PVC, EO	FLL MLS	✓	✓		SLSV025LS		
PVDF Durapel™ Membrane Superhydrophobic for vial venting and transducer protection	0.22	25	N/A	PVC, EO	FLL Vial Spike	✓			SLGVS25LS		
					FLS MLS	✓	✓		SLGVS25XS		
					FLL MLS	✓	✓		SLGVS25PS		
					FLL MLL	✓	✓		SLGVS25US		
MCE MF-Millipore™ Membrane Most referenced general purpose membrane	0.22	25	100 mL (< 100 µL)	PVC, EO	Vented FLL MLS	✓	✓		SLGSV255F		
					FLL MLL	✓	✓		SLGL0250S		
		33	100 mL (< 100 µL)	Modified acrylic, EO	FLL MLS	✓			SLGS033SS		SLGS033SB
						✓	✓		SLGSM33SS		
	0.45	33	100 mL (< 100 µL)	Modified acrylic, EO	FLL MLS	✓			SLHA033SS		SLHA033SB
						✓	✓		SLHAM33SS		
	0.8	25	100 mL (< 100 µL)	PVC, EO	Vented FLL MLS	✓	✓		SLAAV255F		
		33	100 mL (< 100 µL)	Modified acrylic, EO	FLL MLS	✓	✓		SLAA033SS		SLAA033SB
PTFE Fluoropore™ Membrane Hydrophobic for venting and gas filtration	0.2	25	100 mL (< 100 µL)	PVC, EO	FLL MLS	✓	✓		SLFG025LS		
					FLL MLS with Needle	✓		SLFGN25VS			

*Medical devices are only sold in countries where they are registered. Please contact your local customer support team for availability.

Venting and Gas Filtration:

Autoclavable for sterile-filtering gases / Venting sterile containers / In-line vacuum pump protection

Membrane	Pore size (µm)	Diameter (mm)	Housing and Sterilization method	Inlet connection	Outlet connection	10 PK	25 PK	50 PK	100 PK
PTFE Fluoropore™ Membrane Hydrophobic chemistry for gas filtration	0.20 µm	25 mm	PP, Autoclavable	FLL	MLS			SLFG02550	
				HB	HB	SLFG05010			
	0.20 µm	50 mm	PP, Autoclavable	HB (latex)	1/8 in. NPTM	SLFG55010			
					1/8 in. NPTM	SLFG65010			SLFG65000
					1/8 in. NPTM	SLFG75010			SLFG75000
					HB (silicone)	SLFG85010			SLFG85000
	0.45 µm	50 mm	PP, Autoclavable	HB	HB	SLFH05010			SLFH05000
	1.0 µm	50 mm	PP, Autoclavable	HB	HB	SLFA05010			SLFA05000

Non-Sterile Millex® syringe filters you can trust.



Chromatography
HPLC, IC, GC / General particle removal / Industrial / Environmental

Membrane	Diameter (mm)	Pore Size (µm)	Process Volume (hold-up)	50 units/pk	100 units/pk	250 units/pk	1000 units/pk
PTFE Hydrophilic Millipore® LCR Membrane Lowest extractables and excellent solvent resistance	4	0.20	1 mL (< 10 µL)		SLLGR04NL		
		0.45	1 mL (< 10 µL)		SLLHR04NL		
	13	0.20	10 mL (≤15 µL)		SLLGX13NL		SLLGX13NK
		0.45	10 mL (≤15 µL)		SLCRX13NL SLCRX13TL*		SLCRX13NK
	33	0.20	100 mL (≤ 80 µL)	SLLG033NS		SLLG033NB	SLLG033NK
		0.45	100 mL (≤ 80 µL)	SLCR033NS		SLCR033NB	SLCR033NK
PTFE Hydrophilic IC Millex® Filters Low IC extractables/Ion Chromatography Certified	13	0.20	10 mL (≤ 25 µL)		SLLGC13NL		
		0.45	10 mL (≤ 25 µL)		SLLHC13NL		
	25	0.20	100 mL (≤ 100 µL)	SLLGC25NS			
		0.45	100 mL (≤ 100 µL)	SLLHC25NS			
PTFE Hydrophobic Fluoropore™ Membrane Excellent solvent resistance	4	0.20	1 mL (< 10 µL)		SLFGR04NL		
		0.45	1 mL (< 10 µL)		SLFHR04NL		
	13	0.20	10 mL (≤ 15 µL)		SLFGX13NL SLFGX13TL*		SLFGX13NK
		0.45	10 mL (≤ 15 µL)		SLFHX13NL SLFHX13TL*		SLFHX13NK
	33	0.20	100 mL (≤ 80 µL)	SLFG033NS		SLFG033NB	SLFG033NK
		0.45	100 mL (≤ 80 µL)	SLFH033NS		SLFH033NB	SLFH033NK
PVDF Durapore® Membrane Low-protein binding	4	0.22	1 mL (< 10 µL)		SLGVR04NL		SLGVR04NK
		0.45	1 mL (< 10 µL)		SLHVR04NL		SLHVR04NK
	13	0.22	10 mL (≤ 15 µL)		SLGVX13NL SLGVX13TL*		SLGVX13NK
		0.45	10 mL (≤ 15 µL)		SLHVX13NL SLHVX13TL*		SLHVX13NK
	33	0.22	100 mL (≤ 80 µL)	SLGV033NS		SLGV033NB	SLGV033NK
		0.45	100 mL (≤ 80 µL)	SLHV033NS		SLHV033NB	SLHV033NK
Nylon Membrane Broad chemical compatibility	13	0.20	10 mL (≤ 15 µL)		SLGNX13NL SLGNX13TL*		SLGNX13NK
		0.45	10 mL (≤ 15 µL)		SLHNX13NL SLHNX13TL*		SLHNX13NK
	33	0.20	100 mL (≤ 80 µL)	SLGN033NS		SLGN033NB	SLGN033NK
		0.45	100 mL (≤ 80 µL)	SLHN033NS		SLHN033NB	SLHN033NK
PES Millipore Express® Membrane Fastest flow, high throughput	13	0.22	10 mL (≤ 15 µL)		SLGPX13NL		SLGPX13NK
		0.45	10 mL (≤ 15 µL)		SLHPX13NL		SLHPX13NK
	33	0.22	100 mL (≤ 80 µL)	SLGP033NS		SLGP033NB	SLGP033NK
		0.45	100 mL (≤ 80 µL)	SLHP033NS		SLHP033NB	SLHP033NK
MCE General laboratory filtration of aqueous solutions	25	0.22	100 mL (≤ 100 µL)			SLGS025NB†	
		0.45	100 mL (≤ 100 µL)			SLHA025NB†	SLHA025NK†
		0.8	100 mL (≤ 100 µL)			SLAA025NB†	SLAA025NK†

† PVC Housing

Automation Compatible

Dissolution testing / HPLC sample prep

Membrane	Diameter (mm)	Pore Size (µm)	Process Volume (hold-up)	200 (8 x 25) units/pk	1000 units/pk
Durapore® Membrane (PVDF)	25	0.45	100 mL (≤ 100 µL)	SLHVDZ5NZ	SLHVDZ5NK
Millipore® LCR Membrane (Hydrophilic PTFE)	25	0.20	100 mL (≤ 100 µL)	SLLGDZ5NZ	SLLGDZ5NK
		0.45	100 mL (≤ 100 µL)	SLCRDZ5NZ	SLCRDZ5NK
Nylon Membrane	25	0.20	100 mL (≤ 100 µL)	SLGNDZ5NZ	SLGNDZ5NK
		0.45	100 mL (≤ 100 µL)	SLHNDZ5NZ	SLHNDZ5NK
Glass Fiber Filter	25	1.0	100 mL (≤ 250 µL)	SLPBDZ5NZ	SLPBDZ5NK
Durapore® Membrane (PVDF) with glass fiber prefilter	25	0.45	100 mL (≤ 250 µL)	SLHVBZ5NZ	SLHVBZ5NK
Millipore® LCR Membrane (Hydrophilic PTFE) with glass fiber prefilter	25	0.45	100 mL (≤ 250 µL)	SLCRBZ5NZ	SLCRBZ5NK
Nylon Membrane with glass fiber prefilter	25	0.45	100 mL (≤ 250 µL)	SLHNBZ5NZ	SLHNBZ5NK

High Particulate Filtration & Automation Compatible

Chromatography sample preparation / Wine analysis / General particulate removal / Industrial / Environmental

Membrane	Diameter (mm)	Pore Size (µm)	Process Volume (hold-up)	50 units/pk	200 (8 x 25) units/pk	1000 units/pk
Millipore® LCR (Hydrophilic PTFE) with graduated multi-layer glass fiber prefilter	25	0.20	100 mL (< 250 µL)	SLLGM25NS		SLLGM25NK
		0.45	100 mL (< 250 µL)	SLCRM25NS		SLCRM25NK
Durapore® (PVDF) with graduated multi-layer glass fiber prefilter	25	0.45	100 mL (< 250 µL)	SLHVM25NS	SLHVM25NZ	SLHVM25NK
		0.20	100 mL (< 250 µL)	SLGNM25NS		SLGNM25NK
Nylon Membrane with graduated multi-layer glass fiber prefilter	25	0.45	100 mL (< 250 µL)	SLHNM25NS	SLHNM25NZ	SLHNM25NK
		0.20	100 mL (< 250 µL)			

Housing Materials (non-sterile)

Application	4 mm	13 mm	25 mm	33 mm	Millex Simplicity® Filters
Chromatography, automation, high-particulate filtration	High density polyethylene	Polypropylene High density polyethylene	High density polyethylene †PVC where noted	Polypropylene	Polypropylene

Inlet Fittings: Female Luer-Lok™ Outlet Fittings: Male Luer Slip, Male Stepped (4 mm filters), and *Tube outlet where noted.

Simplicity® G2 Filtration System for HPLC Samples (The better way to use Millex® Filters)

A vacuum filtration system that allows simultaneous filtration of multiple samples (1-8) directly into standard sized HPLC, UPLC, GC, sample vials (12 x 32 mm). The Simplicity® G2 filtration system is a faster and easier way to use Millex® 33 mm syringe filters when filtering multiple samples.

Millex® Filters for the Simplicity® G2 System


Membrane	Pore Size (µm)	250 units/pk	1000 units/pk
Hydrophilic PTFE	0.20	SAMP2LGNB	SAMP2LGNK
	0.45	SAMP2LCRB	SAMP2LCRK
Durapore® Membrane (PVDF)	0.22	SAMP2GVNB	SAMP2GVNK
	0.45	SAMP2HVNБ	SAMP2HVNK
Nylon	0.20	SAMP2GNNB	SAMP2GNNK
	0.45	SAMP2HNNB	SAMP2HNNK
PES	0.22	SAMP2GPNB	SAMP2GPNK
	0.45	SAMP2HPNB	SAMP2HPNK

Millex Simplicity® Filters

Membrane	Pore Size (µm)	96 units/pk	384 units/pk
Hydrophilic PTFE	0.20	SAMPLG001	SAMPLG004
	0.45	SAMPLCR01	SAMPLCR04

Simplicity® System and Accessories

Simplicity® G2 System	SAMP2SYSB
Vial Trays	SAMVIALTR
Waste Trays	SAMWASTTR
Tube Set Assembly	SAMTUBING
Replacement Lid	SAMP2LID



The Life Science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, Millex, Durapore, Simplicity, Millex Simplicity and Millipore Express are registered trademarks and Durapel, Fluoropore and MF-Millipore are trademarks of Merck KGaA, Darmstadt, Germany. All other trademarks are the property of their respective owners.

MK_PS2549EN Ver. 5.0 59909 11/2024

