

Certificate of Conformity

HOOD SPP 20G WHITE 35X18CM
113-8202

CHARACTERISTICS:

- Made of non-woven Polypropylene
- White color
- Nature Rubber Latex Free

PHYSICAL PROPERTIES

Material	Non-woven Polypropylene 20grams/m2 white
Sizes	Universal
Break Tensile Properties (N) (MD)	>40
Break Tensile Properties (N) (CD)	>20
Elongation Rate (%) (MD)	>60%
Elongation Rate (%) (CD)	>70%
Height	35cm +/-1cm
Width	38m +/- 1cm
Color	White

PACKING:

100pcs/bag, 5bags/carton

QUALITY ASSURANCE CRITERIA:

Hoods are guaranteed to meet 4.0 AQL level II sampling plan per MIL-std-150E for the following defect criteria:

- | | |
|---------------------------|---------------------------------------|
| 1. Broken elastic band | 4. Holes/tears in fabric material |
| 2. Irregular SPP material | caused by sewing machine |
| 3. Dirty | 5. Physical dimensions non-compliance |
| | 6. Not clean |

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2015**